
















NANO RESEARCH ELEMENTS



CONTENT

GENERAL INFORMATION

→ ABOUT NANO RESEARCH ELEMENTS	02
→ FEATURES	03
→ APPLICATION OF NANOTECHNOLOGY.....	04
→ NANO PARTICLES	05-14
→ MICRO POWDER.....	15-23
→ NANO STRUCTURES.....	24-31
→ QUANTUM DOT'S.....	32-37
→ NANO DISPERSION.....	38-42
→ RARE EARTH MATERIALS	43-44
→ NOVEL MATERIALS	45-46
→ ADVANCE MATERIALS	47-48
→ TARGETS AND WAFERS	49-51
→ CARBON NANOMATERIALS AND GRAPHENE	52

About us,

NanoResearch Elements Inc. is set up by a group of expert Nanotechnologists to concentrate on giving an answer for each mind boggling complex issue with the assistance of Nanotechnology products and applications.

In the most recent decade, improvements in nanotechnology have been fuelled by noteworthy advances in nanomaterials synthesis and in addition, promising application information for the biomedical, electronics, and energy fields. Nanomaterials innovation has the potential to address a number of today's biggest scientific challenges, ranging from the need for more efficient alternative energy technologies, faster and flexible electronics, and improved disease diagnosis and treatment. The distinctive properties of nanomaterials include large surface-to-volume ratios, favourable conductivity, unique optical and confinement effects, as well as many others.

Today, an existence without nanotechnology is difficult to imagine. Nanotechnologies - to be more particular: nanomaterials - are already used in numerous products and industrial applications. Our Nanotechnology Products and Application database as of now give an outline of how nanomaterials and nanostructure applications are utilized today industrial and commercial applications across industries.

NanoResearch Elements offer a huge collection of nanomaterials that navigate the Periodic Table including gold, silver, iron oxide, metal and metal composites, oxides, nitrides and ceramics. We additionally offer an extensive variety of nanowires, nanopowders and nanoparticle dispersions, magnetic nanopowders, nanoparticles and quantum dots. Organic nanomaterials, for example, carbon nanotubes, graphene and graphene oxide, and fullerenes are available with different purity grades and surface modifications. We additionally offer dendrimers, which are organic nanomaterials with well-defined surface chemical groups and sizes suitable for drug delivery research and development.

We trust Nanotechnology will drive real future upgrades for production technology in chemistry, biotechnology, electronics, medicine, material science, alternative energy, lubricants, and agriculture.

FEATURES OF NANO RESEARCH ELEMENTS

**Quality in a service or product is not what you put into it.
It is what the client or customer gets out of it. Our primary goal is client satisfaction.**



INNOVATION

Our interdisciplinary research and development team has a set up record of creating innovative resolution for complex issues.



SCIENTIFIC

We select execution measurements in view of the requirements of the final product, and recognize the basic parameters for the improvement of high return, versatile and reproducible processing and manufacturing techniques.



CONFIDENTIALITY

Our organizational relationship with clients, and in addition all work and information related to the research and development action, will remain completely confidential



ON TIME AND COST EFFECTIVE

Our team is determined to give on time and cost effective deliveries that exceed our customer's expectations.




SINGLE PLATFORM FOR EVERY NEED

We provide Exclusive Variety of Nano products and quality services at single platform



COMPREHENSIVE COMMUNICATION


We formulate documented reports, permit access to lab notebooks, and regular conference calls to answer questions and collaboratively make strategic decisions about the project direction.



CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-1001	ALUMINIUM NANOPARTICLES	Al, 99.9%, APS <100 nm
NRE-1002	ANTIMONY NANOPARTICLES	Sb, 99.9%, APS <100 nm
NRE-1003	BISMUTH NANOPARTICLES	Bi, 99.9%, APS <100 nm
NRE-1004	BORON NANOPARTICLES	B, 99.9%, APS <100 nm
NRE-1005	CARBON NANOPARTICLES	C, 99.9%, APS <100 nm
NRE-1006	CARBONYL IRON NANOPARTICLES	Fe, 99.9%, APS <100 nm
NRE-1007	CHROMIUM NANOPARTICLES	Cr, 99.9%, APS <100 nm
NRE-1008	COBALT NANOPARTICLES	Co, 99.9%, APS <100 nm
NRE-1009	COPPER NANOPARTICLES	Cu, 99.9%, APS <100 nm
NRE-1010	DIAMOND NANOPARTICLES	C, 99.9%, APS <100 nm
NRE-1011	GERMANIUM NANOPARTICLES	Ge, 99.9%, APS <100 nm
NRE-1012	GOLD NANOPARTICLES	Au, 99.9%, APS <100 nm
NRE-1013	HAFNIUM NANOPARTICLES	Hf, 99.9%, APS <100 nm
NRE-1014	INDIUM NANOPARTICLES	In, 99.9%, APS <100 nm
NRE-1015	IRIDIUM NANOPARTICLES	Ir, 99.9%, APS <100 nm
NRE-1016	IRON NANOPARTICLES	Fe, 99.9%, APS <100 nm
NRE-1017	LEAD NANOPARTICLES	Pb, 99.9%, APS <100 nm
NRE-1018	MAGNESIUM NANOPARTICLES	Mg, 99.9%, APS <100 nm
NRE-1019	MANGANESE NANOPARTICLES	Mn, 99.9%, APS <100 nm
NRE-1020	MOLYBDENUM NANOPARTICLES	Mo, 99.9%, APS <100 nm
NRE-1021	NICKEL NANOPARTICLES	Ni, 99.9%, APS <100 nm
NRE-1022	NIObIUM NANOPARTICLES	Nb, 99.9%, APS <100 nm
NRE-1023	OSMIUM NANOPARTICLES	Os, 99.9%, APS <100 nm
NRE-1024	PALLADIUM NANOPARTICLES	Pd, 99.9%, APS <100 nm
NRE-1025	PLATINUM NANOPARTICLES	Pt, 99.9%, APS <100 nm
NRE-1026	RHENIUM NANOPARTICLES	Re, 99.9%, APS <100 nm
NRE-1027	RHODIUM NANOPARTICLES	Rh, 99.9%, APS <100 nm
NRE-1028	RUTHENIUM NANOPARTICLES	Ru, 99.9%, APS <100 nm
NRE-1029	SCANDIUM NANOPARTICLES	Sc, 99.9%, APS <100 nm
NRE-1030	SELENIUM NANOPARTICLES	Se, 99.9%, APS <100 nm
NRE-1031	SILICON NANOPARTICLES	Si, 99.9%, APS <100 nm
NRE-1032	SILVER NANOPARTICLES	Ag, 99.9%, APS <100 nm
NRE-1033	SULPHUR NANOPARTICLES	S, 99.9%, APS <100 nm
NRE-1034	SUPER ACTIVATED CARBON NP'S	C, 99.9%, APS <100 nm
NRE-1035	TANTALUM NANOPARTICLES	Ta, 99.9%, APS <100 nm
NRE-1036	TELLURIUM NANOPARTICLES	Te, 99.9%, APS <100 nm
NRE-1037	TIN NANOPARTICLES	Sn, 99.9%, APS <100 nm
NRE-1038	TITANIUM NANOPARTICLES	Ti, 99.9%, APS <100 nm
NRE-1039	TUNGSTEN NANOPARTICLES	W, 99.9%, APS <100 nm
NRE-1040	VANADIUM NANOPARTICLES	V, 99.9%, APS <100 nm
NRE-1041	YTTRIUM NANOPARTICLES	Y, 99.9%, APS <100 nm
NRE-1042	ZINC NANOPARTICLES	Zn, 99.9%, APS <100 nm
NRE-1043	ZIRCONIUM NANOPARTICLES	Zr, 99.9%, APS <100 nm


CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-2001	ALUMINIUM MAGNESIUM ALLOY NANOPARTICLES	Al-Mg, 99.9%, APS <100 nm
NRE-2002	CHROMIUM CARBIDE NICKEL CHROMIUM ALLOY NP'S	Cr ₃ C ₂ -NiCr, 99.9%, APS <100 nm
NRE-2003	COPPER CADMIUM ALLOY NANOPARTICLES	Cu-Cd, 99.9%, APS <100 nm
NRE-2004	COPPER INDIUM ALLOY NANOPARTICLES	Cu-In, 99.9%, APS <100 nm
NRE-2005	COPPER INDIUM GALLIUM ALLOY NANOPARTICLES	Cu-In-Ga, 99.9%, APS <100 nm
NRE-2006	COPPER INDIUM SULPHUR ALLOY NANOPARTICLES	Cu-In-S, 99.9%, APS <100 nm
NRE-2007	COPPER NICKEL ALLOY NANOPARTICLES	Cu-Ni, 99.9%, APS <100 nm
NRE-2008	COPPER TIN ALLOY NANOPARTICLES	Cu-Sn, 99.9%, APS <100 nm
NRE-2009	COPPER ZINC ALLOY NANOPARTICLES	Cu-Zn, 99.9%, APS <100 nm
NRE-2010	GOLD COPPER ALLOY NANOPARTICLES	Au-Cu, 99.9%, APS <100 nm
NRE-2011	GOLD TIN ALLOY NANOPARTICLES	Au-Sn, 99.9%, APS <100 nm
NRE-2012	IRON BORON ALLOY NANOPARTICLES	Fe-B, 99.9%, APS <100 nm
NRE-2013	IRON CHROMIUM COBALT ALLOY NANOPARTICLES	Fe-Cr-Co, 99.9%, APS <100 nm
NRE-2014	IRON MOLYBDENUM ALLOY NANOPARTICLES	Fe-Mo, 99.9%, APS <100 nm
NRE-2015	IRON NICKEL ALLOY NANOPARTICLES	Fe-Ni, 99.9%, APS <100 nm
NRE-2016	IRON NICKEL COBALT ALLOY NANOPARTICLES	Fe-Ni-Co, 99.9%, APS <100 nm
NRE-2017	IRON NIOBIUM ALLOY NANOPARTICLES	Fe-Nb, 99.9%, APS <100 nm
NRE-2018	IRON SILICON ALLOY NANOPARTICLES	Fe-Si, 99.9%, APS <100 nm
NRE-2019	IRON TITANIUM ALLOY NANOPARTICLES	Fe-Ti, 99.9%, APS <100 nm
NRE-2020	IRON TUNGSTEN ALLOY NANOPARTICLES	Fe-W, 99.9%, APS <100 nm
NRE-2021	IRON VANADIUM ALLOY NANOPARTICLES	Fe-V, 99.9%, APS <100 nm
NRE-2022	LANTHANUM IRON PALLADIUM ALLOY NANOPARTICLES	La-Fe-Pd, 99.9%, APS <100 nm
NRE-2023	NEODYMIUM IRON BORON ALLOY NANOPARTICLES	Nd-Fe-B, 99.9%, APS <100 nm
NRE-2024	NICKEL ALUMINIUM ALLOY NANOPARTICLES	Ni-Al, 99.9%, APS <100 nm
NRE-2025	NICKEL CHROMIUM ALLOY NANOPARTICLES	Ni-Cr, 99.9%, APS <100 nm
NRE-2026	NICKEL CHROMIUM COBALT ALLOY NANOPARTICLES	Ni-Cr-Co, 99.9%, APS <100 nm
NRE-2027	NICKEL GRAPHITE ALLOY NANOPARTICLES	Ni-C, 99.9%, APS <100 nm
NRE-2028	NICKEL TITANIUM ALLOY NANOPARTICLES	Ni-Ti, 99.9%, APS <100 nm
NRE-2029	PALLADIUM SILVER ALLOY NANOPARTICLES	Pd-Ag, 99.9%, APS <100 nm
NRE-2030	PLATINIUM EUROPIUM ALLOY NANOPARTICLES	Pt-Eu, 99.9%, APS <100 nm
NRE-2031	PLATINIUM IRIIDIUM ALLOY NANOPARTICLES	Pt-Ir, 99.9%, APS <100 nm
NRE-2032	SILICON ALUMINIUM ALLOY NANOPARTICLES	Si-Al, 99.9%, APS <100 nm
NRE-2033	SILVER COPPER ALLOY NANOPARTICLES	Ag-Cu, 99.9%, APS <100 nm
NRE-2034	SILVER NICKEL ALLOY NANOPARTICLES	Ag-Ni, 99.9%, APS <100 nm
NRE-2035	SILVER PALLADIUM ALLOY NANOPARTICLES	Ag-Pd, 99.9%, APS <100 nm
NRE-2036	SILVER TIN ALLOY NANOPARTICLES	Ag-Sn, 99.9%, APS <100 nm
NRE-2037	SILVER ZIRCONIUM ALLOY NANOPARTICLES	Ag-Zr, 99.9%, APS <100 nm
NRE-2038	STAINLESS STEEL ALLOY NANOPARTICLES	99.9%, APS <100 nm
NRE-2039	TANTALUM NIOBIUM CARBIDE ALLOY NANOPARTICLES	Ta-NbC, 99.9%, APS <100 nm
NRE-2040	TITANIUM CARBIDE TUNGSTUN CARBIDE ALLOY NP'S	TiC-WC, 99.9%, APS <100 nm
NRE-2041	TITANIUM TANTALUM CARBIDE TUNGSTUN CARBIDE ALLOY NP'S	Ti-TaC-WC, 99.9%, APS <100 nm
NRE-2042	TUNGSTUN CARBIDE CHROMIUM CARBIDE NICKEL ALLOY NP'S	WC-CrC-Ni, 99.9%, APS <100 nm
NRE-2043	TUNGSTUN CARBIDE COBALT ALLOY NANOPARTICLES	WC-Co, 99.9%, APS <100 nm
NRE-2044	TUNGSTUN CARBIDE COBALT CHROMIUM ALLOY NP'S	WC-Co-Cu, 99.9%, APS <100 nm
NRE-2045	TUNGSTUN CARBIDE NICKEL ALLOY NANOPARTICLES	WC-Ni, 99.9%, APS <100 nm

CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-2046	DEVARD'S ALLOY NANOPARTICLES	99.9%, APS <100 nm
NRE-2047	COBALT CHROMIUM ALLOY NANOPARTICLES	99.9%, APS <100 nm
NRE-2048	COBALT CHROMIUM ALUMINIUM YTTRIUM ALLOY NANOPARTICLES	99.9%, APS <100 nm
NRE-2049	NICKEL CHROMIUM ALUMINIUM YTTRIUM ALLOY NANOPARTICLES	99.9%, APS <100 nm
NRE-2050	NICKEL COBALT CHROMIUM ALUMINIUM ALLOY NANOPARTICLES	99.9%, APS <100 nm
NRE-2051	VANADIUM ALUMINIUM ALLOY NANOPARTICLES	99.9%, APS <100 nm
NRE-2052	ZIRCONIUM NICKEL ALLOY NANOPARTICLES	99.9%, APS <100 nm
NRE-2053	ZIRCONIUM ALUMINIUM ALLOY NANOPARTICLES	99.9%, APS <100 nm
NRE-2054	TUNGSTEN CHROMIUM ALLOY NANOPARTICLES	99.9%, APS <100 nm
NRE-2055	TITANIUM ALUMINIDE ALLOY NANOPARTICLES	99.9%, APS <100 nm



CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-3001	ACTIVATED SILICON OXIDE NANOPARTICLES	SiO ₂ , 99.9%, APS <100 nm
NRE-3003	ALUMINIUM OXIDE NANOPARTICLES	Al ₂ O ₃ , 99.9%, APS <100 nm
NRE-3004	ALUMINUM HYDROXIDE NANOPARTICLES	Al(OH) ₃ , 99.9%, APS <100 nm
NRE-3005	ANTIMONY OXIDE NANOPARTICLES	Sb ₂ O ₃ , 99.9%, APS <100 nm
NRE-3006	BISMUTH OXIDE NANOPARTICLES	Bi ₂ O ₃ , 99.9%, APS <100 nm
NRE-3007	CERIUM OXIDE NANOPARTICLES	CeO ₂ , 99.9%, APS <100 nm
NRE-3008	CHROMIUM OXIDE NANOPARTICLES	Cr ₂ O ₃ , 99.9%, APS <100 nm
NRE-3009	COBALT OXIDE NANOPARTICLES	Co ₃ O ₄ , 99.9%, APS 40-80 nm
NRE-3011	COBALT(III) OXIDE NANOPARTICLES	Co ₂ O ₃ , 99.9%, APS <100 nm
NRE-3012	COPPER OXIDE LONG FIBROUS NANOPARTICLES	CuO, 99.9%, APS <100 nm
NRE-3013	COPPER OXIDE NANOPARTICLES	CuO, 99.9%, APS <100 nm
NRE-3014	DYSPROSIUM OXIDE NANOPARTICLES	Dy ₂ O ₃ , 99.9%, APS <100 nm
NRE-3015	EUROPIUM OXIDE NANOPARTICLES	Eu ₂ O ₃ , 99.9%, APS <100 nm
NRE-3016	GADOLINIUM OXIDE NANOPARTICLES	Gd ₂ O ₃ , 99.9%, APS <100 nm
NRE-3017	GALLIUM OXIDE NANOPARTICLES	Ga ₂ O ₃ , 99.9%, APS <100 nm
NRE-3018	HAFNIUM OXIDE NANOPARTICLES	HfO ₂ , 99.9%, APS <100 nm
NRE-3019	INDIUM HYDROXIDE NANOPARTICLES	In(OH) ₃ , 99.9%, APS <100 nm
NRE-3020	INDIUM OXIDE NANOPARTICLES	In ₂ O ₃ , 99.9%, APS <100 nm
NRE-3021	IRON HYDROXIDE NANOPARTICLES	Fe(OH) ₃ , 99.9%, APS <100 nm
NRE-3022	IRON OXIDE NANOPARTICLES	Fe ₂ O ₃ , 99.9%, APS <100 nm
NRE-3023	LANTHANUM OXIDE NANOPARTICLES	La ₂ O ₃ , 99.9%, APS <100 nm
NRE-3025	MAGNESIUM HYDROXIDE NANOPARTICLES	Mg(OH) ₂ , 99.9%, APS <100 nm
NRE-3026	MAGNESIUM OXIDE NANOPARTICLES	MgO, 99.9%, APS <100 nm
NRE-3027	MAGNETITE OXIDE NANOPARTICLES	Fe ₃ O ₄ , 99.9%, APS <100 nm
NRE-3028	MANGANESE DIOXIDE NANOPARTICLES	Mn ₂ O ₃ , 99.9%, APS <100 nm
NRE-3029	MANGANESE OXIDE NANOPARTICLES	Mn ₃ O ₄ , 99.9%, APS <100 nm
NRE-3030	MOLYBDENUM OXIDE NANOPARTICLES	MoO ₃ , 99.9%, APS <100 nm
NRE-3031	NEODYMIUM OXIDE NANOPARTICLES	Nd ₂ O ₃ , 99.9%, APS <100 nm
NRE-3032	NICKEL HYDROXIDE NANOPARTICLES	Ni(OH) ₂ , 99.9%, APS <100 nm
NRE-3033	NICKEL OXIDE NANOPARTICLES	NiO, 99.9%, APS <100 nm
NRE-3034	PRASEODYMIUM OXIDE NANOPARTICLES	Pr ₆ O ₁₁ , 99.9%, APS <100 nm
NRE-3035	SAMARIUM OXIDE NANOPARTICLES	Sm ₂ O ₃ , 99.9%, APS <100 nm
NRE-3036	SILICON OXIDE NANOPARTICLES	SiO ₂ , 99.9%, APS <100 nm
NRE-3037	SILVER OXIDE NANOPARTICLES	Ag ₂ O, 99.9%, APS <100 nm
NRE-3038	SUPERPARAMAGNETIC IRON OXIDE NP'S	Fe ₃ O ₄ , 99.9%, APS <100 nm
NRE-3039	TERBIUM OXIDE NANOPARTICLES	Tb ₄ O ₇ , 99.9%, APS <100 nm
NRE-3040	TIN DIOXIDE NANOPARTICLES	SnO ₂ , 99.9%, APS <100 nm
NRE-3042	TITANIUM OXIDE NANOPARTICLES	TiO ₂ , 99.9%, APS <100 nm
NRE-3043	TUNGSTEN OXIDE NANOPARTICLES	WO ₃ , 99.9%, APS <100 nm
NRE-3044	YTTRIUM OXIDE NANOPARTICLES	Y ₂ O ₃ , 99.9%, APS <100 nm
NRE-3046	ZINC OXIDE NANOPARTICLES	ZnO, 99.9%, APS <100 nm
NRE-3047	ZINC OXIDE NP'S COATED WITH 1WT% SILANE	ZnO, 99.9%, APS <100 nm
NRE-3048	ZIRCONIUM HYDROXIDE NANOPARTICLES	Zr(OH) ₄ , 99.9%, APS <100 nm
NRE-3049	ZIRCONIUM OXIDE NANOPARTICLES	ZrO ₂ , 99.9%, APS <100 nm
NRE-3050	GERMANIUM OXIDE NANOPARTICLES	GeO ₂ , 99.9%, APS <100 nm

CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-4001	ANTIMONY TIN OXIDE (ATO) NANOPARTICLES	$\text{SnO}_2:\text{Sb}_2\text{O}_3$, 99.9%, APS <100 nm
NRE-4002	ALUMINIUM DOPED ZINC (AZO) OXIDE NANOPARTICLES	99.9%, APS <100 nm
NRE-4004	BARIUM IRON OXIDE NANOPARTICLES	$\text{BaFe}_{12}\text{O}_{19}$, 99.9%, APS <100 nm
NRE-4005	BARIUM STRONTIUM TITANATE NANOPARTICLES	BaO_4SrTi , 99.9%, APS <100 nm
NRE-4006	BARIUM TITANATE NANOPARTICLES	BaTiO_3 , 99.9%, APS <100 nm
NRE-4007	COBALT IRON OXIDE NANOPARTICLES	CoFe_2O_4 , 99.9%, APS <100 nm
NRE-4008	INDIUM TIN OXIDE (ITO) NANOPARTICLES	$\text{In}_2\text{O}_3:\text{SnO}_2$, 99.9%, APS <100 nm
NRE-4009	MAGNESIUM ALUMINATE OXIDE NANOPARTICLES	Al_2MgO_4 , 99.9%, APS <100 nm
NRE-4011	MANGANESE IRON OXIDE NANOPARTICLES	MnFe_2O_4 , 99.9%, APS <100 nm
NRE-4012	NICKEL COBALT IRON OXIDE NANOPARTICLES	$\text{Ni}_0.5\text{Co}_0.5\text{Fe}_2\text{O}_4$, 99.9%, APS <100 nm
NRE-4013	NICKEL IRON OXIDE NANOPARTICLES	NiFe_2O_4 , 99.9%, APS <100 nm
NRE-4014	NICKEL ZINC IRON OXIDE NANOPARTICLES	$\text{Ni}_0.5\text{Zn}_0.5\text{Fe}_2\text{O}_4$, 99.9%, APS <100 nm
NRE-4015	PRECIPITATED BARIUM SULFATE OXIDE NANOPARTICLES	BaSO_4 , 99.9%, APS <100 nm
NRE-4016	PRECIPITATED CALCIUM CARBONATE OXIDE NP'S	CaCO_3 , 99.9%, APS <100 nm
NRE-4018	STRONTIUM IRON OXIDE NANOPARTICLES	$\text{SrFe}_{12}\text{O}_{19}$, 99.9%, APS <100 nm
NRE-4019	STRONTIUM TITANATE NANOPARTICLES	SrTiO_3 , 99.9%, APS <100 nm
NRE-4020	YTTRIUM ALUMINATE NANOPARTICLES	$\text{Y}_3\text{Al}_5\text{O}_{12}$, 99.9%, APS <100 nm
NRE-4022	ZINC COBALT IRON OXIDE NANOPARTICLES	$\text{Zn}_0.5\text{Co}_0.5\text{Fe}_2\text{O}_4$, 99.9%, APS <100 nm
NRE-4023	ZINC IRON OXIDE NANOPARTICLES	ZnFe_2O_4 , 99.9%, APS <100 nm
NRE-4024	ZINC MANGANESE IRON OXIDE NANOPARTICLES	$\text{Zn}_0.5\text{Mn}_0.5\text{Fe}_2\text{O}_4$, 99.9%, APS <100 nm
NRE-4025	ZIRCONIA TOUGHENED ALUMINA (ZTA) NANOPARTICLES	$\text{ZrO}_2/\text{Al}_2\text{O}_3$, 99.9%, APS <100 nm
NRE-4026	ZIRCONIA YTTRIA NANOPARTICLES	$\text{ZrO}_2\text{-}3\text{Y}$, 99.9%, APS <100 nm
NRE-4027	ZIRCONIUM BARIUM TITANATE (BZT) NANOPARTICLES	$\text{Zr}_0.2\text{BaTi}_0.8\text{O}_3$, 99.9%, APS <100 nm
NRE-4028	TITANIUM SILICON OXIDE NANOPARTICLES	TiO_2SiO_2 , 99.9%, APS <100 nm




CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-5001	ALUMINIUM NITRIDE NANOPARTICLES	AlN, 99.9%, APS <100 nm
NRE-5002	BARIUM CARBONATE NANOPARTICLES	BaCO ₃ , 99.9%, APS <100 nm
NRE-5003	BARIUM FERRITE NANOPARTICLES	BaFe ₂ O ₇ , 99.9%, APS <100 nm
NRE-5004	BARIUM FLUORIDE NANOPARTICLES	BaF ₂ , 99.9%, APS <100 nm
NRE-5005	BARIUM SULFATE NANOPARTICLES	BaSO ₄ , 99.9%, APS <100 nm
NRE-5006	BARIUM TUNGSTATE NANOPARTICLES	BaWO ₄ , 99.9%, APS <100 nm
NRE-5007	BISMUTH FERRITE NANOPARTICLES	BiFeO ₃ , 99.9%, APS <100 nm
NRE-5008	BISMUTH TELLURIDE NANOPARTICLES	Bi ₂ Te ₃ , 99.9%, APS <100 nm
NRE-5009	BISMUTH TUNGSTATE NANOPARTICLES	Bi ₂ (WO ₄), 99.9%, APS <100 nm
NRE-5010	BORON CARBIDE NANOPARTICLES	B ₄ C, 99.9%, APS <100 nm
NRE-5011	BORON NITRIDE NANOPARTICLES	BN, 99.9%, APS <100 nm
NRE-5012	BORON SILICIDE NANOPARTICLES	SiB ₆ , 99.9%, APS <100 nm
NRE-5013	CADMIUM NITRATE NANOPARTICLES	Cd(NO ₃) ₂ , 99.9%, APS <100 nm
NRE-5014	CADMIUM SULFIDE NANOPARTICLES	CdS, 99.9%, APS <100 nm
NRE-5015	CADMIUM TELLURIDE NANOPARTICLES	CdTe, 99.9%, APS <100 nm
NRE-5016	CADMIUM TUNGSTATE NANOPARTICLES	CdWO ₄ , 99.9%, APS <100 nm
NRE-5017	CALCIUM CARBONATE NANOPARTICLES	CaCO ₃ , 99.9%, APS <100 nm
NRE-5018	CALCIUM SULFATE NANOPARTICLES	CaSO ₄ , 99.9%, APS <100 nm
NRE-5019	CALCIUM SULFIDE NANOPARTICLES	CaF ₂ , 99.9%, APS <100 nm
NRE-5020	CARBON ALUMINUM NITRIDE NANOPARTICLES	AlNC, 99.9%, APS <100 nm
NRE-5021	CARBON TITANIUM NITRIDE NANOPARTICLES	TiNC, 99.9%, APS <100 nm
NRE-5022	CESIUM CARBONATE NANOPARTICLES	Ce ₂ CO ₃ , 99.9%, APS <100 nm
NRE-5023	CESIUM TELLURIDE NANOPARTICLES	Ce ₂ TeO ₃ , 99.9%, APS <100 nm
NRE-5024	CESIUM TUNGSTEN NANOPARTICLES	Ce ₂ WO ₄ , 99.9%, APS <100 nm
NRE-5025	CHROMIUM CARBIDE NANOPARTICLES	Cr ₃ C ₂ , 99.9%, APS <100 nm
NRE-5026	CHROMIUM NITRIDE NANOPARTICLES	CrN, 99.9%, APS <100 nm
NRE-5027	COBALT SULPHIDE NANOPARTICLES	Co _x S _y , 99.9%, APS <100 nm
NRE-5028	COBALT ZINC FERRITE NANOPARTICLES	Co _{0.5} Zn _{0.5} Fe ₂ O ₄ , 99.9%, APS <100 nm
NRE-5029	COPPER MONOSULFIDE NANOPARTICLES	CuS, 99.9%, APS <100 nm
NRE-5030	COPPER SILICATE NANOPARTICLES	CuO ₃ Si, 99.9%, APS <100 nm
NRE-5031	GALLIUM SULFATE NANOPARTICLES	Ga ₂ (SO ₄), 99.9%, APS <100 nm
NRE-5032	GALLIUM TRICHLORIDE NANOPARTICLES	GaCl ₃ , 99.9%, APS <100 nm
NRE-5033	GOLD CHLORIDE NANOPARTICLES	AuCl, 99.9%, APS <100 nm
NRE-5034	HAFNIUM CARBIDE NANOPARTICLES	HfC, 99.9%, APS <100 nm
NRE-5035	HAFNIUM CHLORIDE NANOPARTICLES	HfCl ₄ , 99.9%, APS <100 nm
NRE-5036	HAFNIUM DIBORIDE NANOPARTICLES	HfB ₂ , 99.9%, APS <100 nm
NRE-5037	HYDRATED LIME NANOPARTICLES	Ca(OH) ₂ , 99.9%, APS <100 nm
NRE-5038	INDIUM CHLORIDE NANOPARTICLES	InCl ₃ , 99.9%, APS <100 nm
NRE-5039	INDIUM FLUORIDE NANOPARTICLES	InF ₃ , 99.9%, APS <100 nm
NRE-5040	INDIUM NITRIDE NANOPARTICLES	In ₃ N ₃ ·xH ₂ O, 99.9%, APS <100 nm
NRE-5041	INDIUM SULFATE NANOPARTICLES	In ₂ (SO ₄) ₃ , 99.9%, APS <100 nm
NRE-5042	INDIUM SULFIDE NANOPARTICLES	In ₂ S ₃ , 99.9%, APS <100 nm
NRE-5043	INDIUM TELLURIDE NANOPARTICLES	In ₂ Te ₃ , 99.9%, APS <100 nm
NRE-5044	IRIDIUM CHLORIDE NANOPARTICLES	IrCl ₃ , 99.9%, APS <100 nm
NRE-5045	IRON BROMIDE NANOPARTICLES	FeBr ₂ , 99.9%, APS <100 nm

CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-5046	IRON CARBIDE NANOPARTICLES	Fe ₅ C ₂ , 99.9%, APS <100 nm
NRE-5047	IRON SULFIDE NANOPARTICLES	FeS, 99.9%, APS <100 nm
NRE-5048	LANTHANUM CHLORIDE NANOPARTICLES	LaCl ₃ , 99.9%, APS <100 nm
NRE-5049	LANTHANUM FLUORIDE NANOPARTICLES	LaF ₃ , 99.9%, APS <100 nm
NRE-5050	LANTHANUM HEXABORIDE NANOPARTICLES	LaB ₆ , 99.9%, APS <100 nm
NRE-5051	LANTHANUM TRIFLUORIDE NANOPARTICLES	LaF ₃ , 99.9%, APS <100 nm
NRE-5052	LEAD NITRIDE NANOPARTICLES	Pb(NO ₃) ₂ , 99.9%, APS <100 nm
NRE-5053	LEAD SULFIDE NANOPARTICLES	PbS, 99.9%, APS <100 nm
NRE-5054	LEAD TELLURIDE NANOPARTICLES	PbTe, 99.9%, APS <100 nm
NRE-5055	LEAD TUNGSTEN NANOPARTICLES	PbWO ₄ , 99.9%, APS <100 nm
NRE-5056	LITHIUM BOROHYDRIDE NANOPARTICLES	LiBH ₄ , 99.9%, APS <100 nm
NRE-5057	LITHIUM FLUORIDE NANOPARTICLES	LiF, 99.9%, APS <100 nm
NRE-5058	LITHIUM TELLURIDE NANOPARTICLES	Li ₂ Te, 99.9%, APS <100 nm
NRE-5059	LITHIUM TUNGSTEN NANOPARTICLES	Li ₂ WO ₄ , 99.9%, APS <100 nm
NRE-5060	MAGNESIUM CARBONATE NANOPARTICLES	MgCO ₃ , 99.9%, APS <100 nm
NRE-5061	MAGNESIUM CHLORIDE NANOPARTICLES	MgCl ₂ , 99.9%, APS <100 nm
NRE-5062	MAGNESIUM FLURIDE NANOPARTICLES	MgF ₂ , 99.9%, APS <100 nm
NRE-5063	MAGNESIUM NITRATE NANOPARTICLES	Mg(NO ₃) ₂ , 99.9%, APS <100 nm
NRE-5064	MAGNESIUM NITRIDE NANOPARTICLES	Mg ₃ N ₂ , 99.9%, APS <100 nm
NRE-5065	MAGNESIUM TUNGSTEN NANOPARTICLES	MgWO ₄ , 99.9%, APS <100 nm
NRE-5066	MANGANESE NITRIDE NANOPARTICLES	Mn ₃ N ₂ , 99.9%, APS <100 nm
NRE-5067	MANGANESE SULPHIDE NANOPARTICLES	MnS, 99.9%, APS <100 nm
NRE-5068	MANGANESE TUNGSTEN NANOPARTICLES	MnO ₄ W, 99.9%, APS <100 nm
NRE-5069	MOLYBDENUM CARBIDE NANOPARTICLES	Mo ₂ C, 99.9%, APS <100 nm
NRE-5070	MOLYBDENUM DISILICIDE NANOPARTICLES	MoSi ₂ , 99.9%, APS <100 nm
NRE-5071	MOLYBDENUM DISULFIDE NANOPARTICLES	MoS ₂ , 99.9%, APS <100 nm
NRE-5072	NICKEL ZINC FERRITE NANOPARTICLES	Ni _{0.5} Zn _{0.5} Fe ₂ O ₄ , 99.9%, APS <100 nm
NRE-5073	NIوبيUM CARBIDE NANOPARTICLES	NbC, 99.9%, APS <100 nm
NRE-5074	POTASSIUM SULFATE NANOPARTICLES	K ₂ SO ₄ , 99.9%, APS <100 nm
NRE-5075	POTASSIUM TELLURITE NANOPARTICLES	K ₂ Te, 99.9%, APS <100 nm
NRE-5076	POTASSIUM TITANATE NANOPARTICLES	K ₂ O.8TiO ₂ , 99.9%, APS <100 nm
NRE-5077	POTASSIUM TUNGSTATE NANOPARTICLES	K ₂ WO ₄ , 99.9%, APS <100 nm
NRE-5078	PTFE NANOPARTICLES	(C ₂ F ₄) _n , 99.9%, APS <100 nm
NRE-5079	RUBIDIUM TUNGSTATE NANOPARTICLES	Rb ₂ WO ₄ , 99.9%, APS <100 nm
NRE-5080	Si ₃ N ₄ NANOPARTICLES	Si ₃ N ₄ , 99.9%, APS <100 nm
NRE-5081	SILICON CARBIDE NANOPARTICLES	SiC, 99.9%, APS <100 nm
NRE-5082	SILICON NITRIDE NANOPARTICLES	Si ₃ N ₄ , 99.9%, APS <100 nm
NRE-5083	SILVER FLUORIDE NANOPARTICLES	AgF, 99.9%, APS <100 nm
NRE-5084	SILVER TUNGSTATE NANOPARTICLES	Ag ₂ WO ₄ , 99.9%, APS <100 nm
NRE-5085	SODIUM FLURIDE NANOPARTICLES	NaF, 99.9%, APS <100 nm
NRE-5086	STRONTIUM CARBONATE NANOPARTICLES	SrCO ₃ , 99.9%, APS <100 nm
NRE-5087	STRONTIUM FLUORIDE NANOPARTICLES	SrF ₂ , 99.9%, APS <100 nm
NRE-5088	STRONTIUM TITANATE NANOPARTICLES	SrTiO ₃ , 99.9%, APS <100 nm
NRE-5089	STRONTIUM TUNGSTATE NANOPARTICLES	SrWO ₄ , 99.9%, APS <100 nm
NRE-5090	TANTALUM CARBIDE NANOPARTICLES	TaC, 99.9%, APS <100 nm

CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-5091	TANTALUM NIOBIUM CARBIDE NP'S	TaC:NbC=90:10, 99.9%, APS <100 nm
NRE-5092	TELLURIUM CHLORIDE NANOPARTICLES	Cl ₄ Te, 99.9%, APS <100 nm
NRE-5093	THALLIUM CARBONATE NANOPARTICLES	Tl ₂ CO ₃ , 99.9%, APS <100 nm
NRE-5094	THALLIUM NITRATE NANOPARTICLES	TlNO ₃ , 99.9%, APS <100 nm
NRE-5095	THALLIUM SULFATE NANOPARTICLES	Tl ₂ SO ₄ , 99.9%, APS <100 nm
NRE-5096	TIN CHLORIDE NANOPARTICLES	SnCl ₄ , 99.9%, APS <100 nm
NRE-5097	TIN II SULFIDE NANOPARTICLES	SnS, 99.9%, APS <100 nm
NRE-5098	TIN NITRIDE NANOPARTICLES	Sn (NO ₃) ₄ , 99.9%, APS <100 nm
NRE-5099	TIN SELENIDE NANOPARTICLES	SnSe, 99.9%, APS <100 nm
NRE-5100	TITANIUM BORIDE NANOPARTICLES	TiB ₂ , 99.9%, APS <100 nm
NRE-5101	TITANIUM CARBIDE NANOPARTICLES	TiC, 99.9%, APS <100 nm
NRE-5102	TITANIUM DIBORIDE NANOPARTICLES	TiB ₂ , 99.9%, APS <100 nm
NRE-5103	TITANIUM DISILICIDE NANOPARTICLES	TiSi ₂ , 99.9%, APS <100 nm
NRE-5104	TITANIUM HYDRIDE NANOPARTICLES	TiH ₂ , 99.9%, APS <100 nm
NRE-5105	TITANIUM NITRIDE NANOPARTICLES	TiSiO ₄ , 99.9%, APS <100 nm
NRE-5107	TUNGSTEN CARBIDE COBALT NANOPARTICLES	WCCo, 99.9%, APS <100 nm
NRE-5108	TUNGSTEN CARBIDE NANOPARTICLES	WC, 99.9%, APS <100 nm
NRE-5109	TUNGSTEN DISULFIDE NANOPARTICLES	WS ₂ , 99.9%, APS <100 nm
NRE-5110	TUNGSTEN TITANIUM CARBIDE NANOPARTICLES	WC:TiC=50:50, 99.9%, APS <100 nm
NRE-5111	TUNGSTEN TITANIUM TANTALUM CARBIDE NP'S	WC:TiC:TaC=50:30:20, 99.9%, APS <100 nm
NRE-5112	VANADIUM CARBIDE NANOPARTICLES	VC, 99.9%, APS <100 nm
NRE-5113	VANADIUM SULFATE NANOPARTICLES	V ₂ (SO ₄) ₃ , 99.9%, APS <100 nm
NRE-5114	ZINC CARBONATE NANOPARTICLES	ZnCO ₃ , 99.9%, APS <100 nm
NRE-5115	ZINC CHLORIDE NANOPARTICLES	ZnCl ₂ , 99.9%, APS <100 nm
NRE-5118	ZINC SULPHIDE NANOPARTICLES	ZnS, 99.9%, APS <100 nm
NRE-5119	ZINC TITANATE NANOPARTICLES	ZnTiO ₃ , 99.9%, APS <100 nm
NRE-5120	ZIRCONIUM (II) HYDRIDE NANOPARTICLES	ZrH ₂ , 99.9%, APS <100 nm
NRE-5121	ZIRCONIUM CARBIDE NANOPARTICLES	ZrC, 99.9%, APS <100 nm
NRE-5122	ZIRCONIUM DIBORIDE NANOPARTICLES	ZrB ₂ , 99.9%, APS <100 nm
NRE-5123	ZIRCONIUM DISILICIDE NANOPARTICLES	ZrSi ₂ , 99.9%, APS <100 nm
NRE-5124	ZIRCONIUM NITRIDE NANOPARTICLES	ZrN, 99.9%, APS <100 nm
NRE-5125	ZIRCONIUM TUNGSTATE NANOPARTICLES	Zr(WO ₄) ₂ , 99.9%, APS <100 nm
NRE-5126	VANADIUM HYDRIDE NANOPARTICLES	H ₅ V-5, 99.9%, APS <100 nm
NRE-5127	HAFNIUM HYDRIDE NANOPARTICLES	HfH ₂ , 99.9%, APS <100 nm
NRE-5128	CALCIUM HYDRIDE NANOPARTICLES	CaH ₂ , 99.9%, APS <100 nm
NRE-5129	ALUMINIUM CARBIDE NANOPARTICLES	Al ₄ C ₃ , 99.9%, APS <100 nm
NRE-5130	MANGANESE CARBIDE NANOPARTICLES	C ₂ Mg, 99.9%, APS <100 nm
NRE-5131	TITANIUM CARBONITRIDE NANOPARTICLES	TiCN, 99.9%, APS <100 nm
NRE-5132	CHROMIUM BORIDE NANOPARTICLES	Cr ₃ B ₄ , 99.9%, APS <100 nm
NRE-5133	NICKEL BORIDE NANOPARTICLES	Ni ₂ B, 99.9%, APS <100 nm
NRE-5134	MAGNESIUM BORIDE NANOPARTICLES	MgB ₂ , 99.9%, APS <100 nm
NRE-5135	SILICON BORIDE NANOPARTICLES	SiB ₆ , 99.9%, APS <100 nm
NRE-5136	ALUMINIUM BORIDE NANOPARTICLES	AlB ₂ , 99.9%, APS <100 nm
NRE-5137	MOLYBDENUM BORIDE NANOPARTICLES	BMo, 99.9%, APS <100 nm
NRE-5138	MANGANESE BORIDE NANOPARTICLES	MgB ₂ , 99.9%, APS <100 nm

CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-5139	NIOBIUM BORIDE NANOPARTICLES	NbB ₂ , 99.9%, APS <100 nm
NRE-5140	COBALT BORIDE NANOPARTICLES	Co ₂ B, 99.9%, APS <100 nm
NRE-5141	VANADIUM BORIDE NANOPARTICLES	B ₂ V, 99.9%, APS <100 nm
NRE-5142	TUNGSTEN BORIDE NANOPARTICLES	B ₅ W ₂ , 99.9%, APS <100 nm
NRE-5143	CALCIUM NITRIDE NANOPARTICLES	Ca ₃ N ₂ , 99.9%, APS <100 nm
NRE-5144	HAFNIUM NITRIDE NANOPARTICLES	HfN, 99.9%, APS <100 nm
NRE-5145	ALUMINIUM SILICIDE NANOPARTICLES	Al ₄ Si ₃ , 99.9%, APS <100 nm
NRE-5146	CHROMIUM SILICIDE NANOPARTICLES	CrSi ₂ , 99.9%, APS <100 nm
NRE-5147	MAGNESIUM SILICIDE NANOPARTICLES	Mg ₂ Si, 99.9%, APS <100 nm
NRE-5148	CALCIUM SILICIDE NANOPARTICLES	CaSi ₂ , 99.9%, APS <100 nm
NRE-5149	VANADIUM SILICIDE NANOPARTICLES	VS ₂ , 99.9%, APS <100 nm
NRE-5150	MANGANESE SILICIDE NANOPARTICLES	MnSi ₂ , 99.9%, APS <100 nm
NRE-5151	FERRIC SILICIDE NANOPARTICLES	FeSi, 99.9%, APS <100 nm
NRE-5152	NIOBIUM SILICIDE NANOPARTICLES	NbSi ₂ , 99.9%, APS <100 nm
NRE-5153	TUNGSTEN SILICIDE NANOPARTICLES	WSi ₂ , 99.9%, APS <100 nm
NRE-5154	NICKEL SILICIDE NANOPARTICLES	Ni ₂ Si, 99.9%, APS <100 nm



CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-6001	CELLULOSE NANOPARTICLES	99%+, APS <100 nm
NRE-6002	CHITOSAN NANOPARTICLES	99%+, APS <100 nm
NRE-6003	CLAY NANOPARTICLES	99%+, APS <100 nm
NRE-6004	EXPANDED PERLITE NANOPARTICLES	99%+, APS <100 nm
NRE-6005	GELATIN NANOPARTICLES	99%+, APS <100 nm
NRE-6006	HALLOYSITE CLAY NANOTUBES	99%+, APS <100 nm
NRE-6007	HYDROXYAPATITE NANOPARTICLES	99%+, APS <100 nm
NRE-6008	KAOLIN NANOCLAY	99%+, APS <100 nm
NRE-6009	MICA NANOPARTICLES	99%+, APS <100 nm
NRE-6010	MONTMORILLONITE BENTONITE NANOCLAY	99%+, APS <100 nm
NRE-6011	MONTMORILLONITE NANOCLAY	99%+, APS <100 nm
NRE-6012	NANOCLAY MODIFIED ASPHALT MATERIALS	99%+, APS <100 nm
NRE-6013	NANOCLAY POWDER	99%+, APS <100 nm
NRE-6014	PERLITE CLAY MIX	99%+, APS <100 nm
NRE-6015	PERLITE NANOCLAY POWDER	99%+, APS <100 nm
NRE-6016	SAPONITE NANOCLAY	99%+, APS <100 nm
NRE-6017	SILVER COATED HYDROXYAPATITE NP'S	99%+, APS <100 nm
NRE-6018	ZEOLITE NANOPARTICLES	99%+, APS <100 nm
NRE-6019	SUPER GRADE MICA POWDER	99%+, APS <100 nm
NRE-6020	CALCINED MICA POWDER	99%+, APS <100 nm
NRE-6021	BENTONITE NANOCLAY	99%+, APS <100 nm



CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-8001	ALUMINUM MICRO POWDER	Al, 99.9%, APS < 40 µm
NRE-8002	ANTIMONY MICRO POWDER	Sb, 99.9%, APS < 40 µm
NRE-8003	BISMUTH MICRO POWDER	Bi, 99.9%, APS < 40 µm
NRE-8004	BORON MICRO POWDER	B, 99.9%, APS < 40 µm
NRE-8005	CARBON MICRO POWDER	C, 99.9%, APS < 40 µm
NRE-8006	CHROMIUM MICRO POWDER	Cr, 99.9%, APS < 40 µm
NRE-8007	COBALT MICRO POWDER	Co, 99.9%, APS < 40 µm
NRE-8008	COPPER MICRO POWDER	Cu, 99.9%, APS < 40 µm
NRE-8009	GERMANIUM MICRO POWDER	Ge, 99.9%, APS < 40 µm
NRE-8010	GOLD MICRO POWDER	Au, 99.9%, APS < 40 µm
NRE-8011	GRAPHITE MICRO POWDER	C, 99.9%, APS < 40 µm
NRE-8012	HAFNIUM MICRO POWDER	Hf, 99.9%, APS < 40 µm
NRE-8013	INDIUM MICRO POWDER	In, 99.9%, APS < 40 µm
NRE-8014	IRIDIUM MICRO POWDER	Ir, 99.9%, APS < 40 µm
NRE-8015	IRON MICRO POWDER	Fe, 99.9%, APS < 40 µm
NRE-8016	CARBONYL IRON MICRO POWDER	Fe, 99.9%, APS < 40 µm
NRE-8017	LEAD MICRO POWDER	Pb, 99.9%, APS < 40 µm
NRE-8018	MAGNESIUM MICRO POWDER	Mg, 99.9%, APS < 40 µm
NRE-8019	MANGANESE MICRO POWDER	Mn, 99.9%, APS < 40 µm
NRE-8020	MOLYBDENUM MICRO POWDER	Mo, 99.9%, APS < 40 µm
NRE-8021	NICKEL MICRO POWDER	Ni, 99.9%, APS < 40 µm
NRE-8022	NIOBIUM MICRO POWDER	Nb, 99.9%, APS < 40 µm
NRE-8023	OSMIUM MICRO POWDER	Os, 99.9%, APS < 40 µm
NRE-8024	PALLADIUM MICRO POWDER	Pd, 99.9%, APS < 40 µm
NRE-8025	PLATINUM MICRO POWDER	Pt, 99.9%, APS < 40 µm
NRE-8026	RHENIUM MICRO POWDER	Re, 99.9%, APS < 40 µm
NRE-8027	RHODIUM MICRO POWDER	Rh, 99.9%, APS < 40 µm
NRE-8028	RUTHENIUM MICRO POWDER	Ru, 99.9%, APS < 40 µm
NRE-8029	SCANDIUM MICRO POWDER	Sc, 99.9%, APS < 40 µm
NRE-8030	SELENIUM MICRO POWDER	Se, 99.9%, APS < 40 µm
NRE-8031	SILICON MICRO POWDER	Si, 99.9%, APS < 40 µm
NRE-8032	SILVER MICRO POWDER	Ag, 99.9%, APS < 40 µm
NRE-8033	SULFUR MICRO POWDER	S, 99.9%, APS < 40 µm
NRE-8034	TANTALUM MICRO POWDER	Ta, 99.9%, APS < 40 µm
NRE-8035	TELLURIUM MICRO POWDER	Te, 99.9%, APS < 40 µm
NRE-8036	TIN MICRO POWDER	Sn, 99.9%, APS < 40 µm
NRE-8037	TITANIUM MICRO POWDER	Ti, 99.9%, APS < 40 µm
NRE-8038	TUNGSTEN MICRO POWDER	W, 99.9%, APS < 40 µm
NRE-8039	VANADIUM MICRO POWDER	V, 99.9%, APS < 40 µm
NRE-8040	YTTRIUM MICRO POWDER	Y, 99.9%, APS < 40 µm
NRE-8041	ZINC MICRO POWDER	Zn, 99.9%, APS < 40 µm
NRE-8042	ZIRCONIUM MICRO POWDER	Zr, 99.9%, APS < 40 µm

CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-9001	ALUMINIUM MAGNESIUM ALLOY MICRO POWDER	Al-Mg, 99.9%, APS <40 µm
NRE-9002	CHROMIUM CARBIDE NICKEL CHROMIUM ALLOY POWDER	Cr3C2-NiCr, 99.9%, APS <40 µm
NRE-9003	COPPER CADMIUM ALLOY MICRO POWDER	Cu-Cd, 99.9%, APS <40 µm
NRE-9004	COPPER INDIUM ALLOY MICRO POWDER	Cu-In, 99.9%, APS <40 µm
NRE-9005	COPPER INDIUM GALLIUM ALLOY MICRO POWDER	Cu-In-Ga, 99.9%, APS <40 µm
NRE-9006	COPPER INDIUM SULPHUR ALLOY MICRO POWDER	Cu-In-S, 99.9%, APS <40 µm
NRE-9007	COPPER NICKEL ALLOY MICRO POWDER	Cu-Ni, 99.9%, APS <40 µm
NRE-9008	COPPER TIN ALLOY MICRO POWDER	Cu-Sn, 99.9%, APS <40 µm
NRE-9009	COPPER ZINC ALLOY MICRO POWDER	Cu-Zn, 99.9%, APS <40 µm
NRE-9010	COPPER ZIRCONIUM ALLOY MICRO POWDER	Cu-Zr, 99.9%, APS <40 µm
NRE-9011	GOLD COPPER ALLOY MICRO POWDER	Au-Cu, 99.9%, APS <40 µm
NRE-9012	GOLD TIN ALLOY MICRO POWDER	Au-Sn, 99.9%, APS <40 µm
NRE-9013	IRON BORON ALLOY MICRO POWDER	Fe-B, 99.9%, APS <40 µm
NRE-9014	IRON CHROMIUM COBALT ALLOY MICRO POWDER	Fe-Cr-Co, 99.9%, APS <40 µm
NRE-9015	IRON MOLYBDENUM ALLOY MICRO POWDER	Fe-Mo, 99.9%, APS <40 µm
NRE-9016	IRON NICKEL ALLOY MICRO POWDER	Fe-Ni, 99.9%, APS <40 µm
NRE-9017	IRON NICKEL COBALT ALLOY MICRO POWDER	Fe-Ni-Co, 99.9%, APS <40 µm
NRE-9018	IRON NIOBIUM ALLOY MICRO POWDER	Fe-Nb, 99.9%, APS <40 µm
NRE-9019	IRON SILICON ALLOY MICRO POWDER	Fe-Si, 99.9%, APS <40 µm
NRE-9020	IRON TITANIUM ALLOY MICRO POWDER	Fe-Ti, 99.9%, APS <40 µm
NRE-9021	IRON TUNGSTEN ALLOY MICRO POWDER	Fe-W, 99.9%, APS <40 µm
NRE-9022	IRON VANADIUM ALLOY MICRO POWDER	Fe-V, 99.9%, APS <40 µm
NRE-9023	LANTHANUM IRON PALLADIUM ALLOY MICRO POWDER	La-Fe-Pd, 99.9%, APS <40 µm
NRE-9024	NEODYMIUM IRON BORON ALLOY MICRO POWDER	Nd-Fe-B, 99.9%, APS <40 µm
NRE-9025	NICKEL ALUMINIUM ALLOY MICRO POWDER	Ni-Al, 99.9%, APS <40 µm
NRE-9026	NICKEL CHROMIUM ALLOY MICRO POWDER	Ni-Cr, 99.9%, APS <40 µm
NRE-9028	NICKEL CHROMIUM COBALT ALLOY MICRO POWDER	Ni-Cr-Co, 99.9%, APS <40 µm
NRE-9029	NICKEL GRAPHITE ALLOY MICRO POWDER	Ni-C, 99.9%, APS <40 µm
NRE-9030	NICKEL TITANIUM ALLOY MICRO POWDER	Ni-Ti, 99.9%, APS <40 µm
NRE-9031	PALLADIUM SILVER ALLOY MICRO POWDER	Pd-Ag, 99.9%, APS <40 µm
NRE-9032	PLATINIUM EUROPIUM ALLOY MICRO POWDER	Pt-Eu, 99.9%, APS <40 µm
NRE-9033	PLATINIUM IRIIDIUM ALLOY MICRO POWDER	Pt-Ir, 99.9%, APS <40 µm
NRE-9034	SILICON ALUMINIUM ALLOY MICRO POWDER	Si-Al, 99.9%, APS <40 µm
NRE-9035	SILVER COPPER ALLOY MICRO POWDER	Ag-Cu, 99.9%, APS <40 µm
NRE-9036	SILVER NICKEL ALLOY MICRO POWDER	Ag-Ni, 99.9%, APS <40 µm
NRE-9037	SILVER PALLADIUM ALLOY MICRO POWDER	Ag-Pd, 99.9%, APS <40 µm
NRE-9038	SILVER TIN ALLOY MICRO POWDER	Ag-Sn, 99.9%, APS <40 µm
NRE-9039	SILVER ZIRCONIUM ALLOY MICRO POWDER	Ag-Zr, 99.9%, APS <40 µm
NRE-9040	STAINLESS STEEL ALLOY MICRO POWDER	99.9%, APS <40 µm
NRE-9041	TITANIUM NIOBIUM CARBIDE ALLOY MICRO POWDER	Ti-NbC, 99.9%, APS <40 µm
NRE-9042	TITANIUM CARBIDE TUNGSTEN CARBIDE ALLOY POWDER	TiC-WC, 99.9%, APS <40 µm
NRE-9043	TITANIUM TANTALUM CARBIDE TUNGSTEN CARBIDE ALLOY	Ti-TaC-WC, 99.9%, APS <40 µm
NRE-9044	TUNGSTUN CARBIDE CHROMIUM CARBIDE NICKEL ALLOY	WC-CrC-Ni, 99.9%, APS <40 µm
NRE-9045	TUNGSTUN CARBIDE COBALT ALLOY MICRO POWDER	WC-CO, 99.9%, APS <40 µm
NRE-9046	TUNGSTUN CARBIDE COBALT CHROMIUM ALLOY MICRO POWDER	WC-Co-Cr, 99.9%, APS <40 µm

Micro Powder

Alloy Micro Powder

CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-9047	TUNGSTUN CARBIDE NICKEL ALLOY MICRO POWDER	99.9%, APS <40 µm
NRE-9048	STAINLESS STEEL 304L MICRO POWDER	99.9%, APS <40 µm
NRE-9050	STAINLESS STEEL 316L MICRO POWDER	99.9%, APS <40 µm
NRE-9051	STAINLESS STEEL 430L MICRO POWDER	99.9%, APS <40 µm
NRE-9052	VANADIUM ALUMINIUM ALLOY POWDER	99.9%, APS<40µm
NRE-9053	ZIRCONIUM NICKEL ALLOY POWDER	99.9%, APS<40µm
NRE-9054	ZIRCONIUM ALUMINIUM ALLOY POWDER	99.9%, APS<40µm
NRE-9055	TUNGSTEN CHROMIUM ALLOY POWDER	99.9%, APS<40µm
NRE-9056	TITANIUM ALUMINIDE ALLOY POWDER	TiAl3, 99.9%, APS<40µm



[illegible]

CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-11001	AZO MICRO POWDER	99.9%, APS<40μm
NRE-11002	BARIUM FLUORIDE MICRO POWDER	BaF ₂ , 99.9%, APS<40μm
NRE-11003	BARIUM SULFATE MICRO POWDER	BaSO ₄ , 99.9%, APS<40μm
NRE-11004	BARIUM TUNGSTATE MICRO POWDER	BaWO ₄ , 99.9%, APS<40μm
NRE-11005	BISMUTH FERRITE MICRO POWDER	BiFeO ₃ , 99.9%, APS<40μm
NRE-11006	BISMUTH TELLURIDE MICRO POWDER	Bi ₂ Te ₃ , 99.9%, APS<40μm
NRE-11007	BISMUTH TUNGSTATE MICRO POWDER	Bi ₂ (WO ₄), 99.9%, APS<40μm
NRE-11008	BORON CARBIDE MICRO POWDER	B ₄ C, 99.9%, APS<40μm
NRE-11009	BORON NITRIDE MICRO POWDER	BN, 99.9%, APS<40μm
NRE-11010	BORON SILICIDE MICRO POWDER	SiB ₆ , 99.9%, APS<40μm
NRE-11011	CADMIUM NITRATE MICRO POWDER	Cd(NO ₃) ₂ , 99.9%, APS<40μm
NRE-11012	CADMIUM SULFIDE MICRO POWDER	CdS, 99.9%, APS<40μm
NRE-11013	CADMIUM TELLURIDE MICRO POWDER	CdTe, 99.9%, APS<40μm
NRE-11014	CADMIUM TUNGSTATE MICRO POWDER	CdWO ₄ , 99.9%, APS<40μm
NRE-11015	CALCIUM BORIDE MICRO POWDER	CaB ₆ , 99.9%, APS<40μm
NRE-11016	CALCIUM SULFIDE MICRO POWDER	CaS, 99.9%, APS<40μm
NRE-11017	CARBON ALUMINUM NITRIDE MICRO POWDER	AlNC, 99.9%, APS<40μm
NRE-11018	CARBON TITANIUM NITRIDE MICRO POWDER	TiNC, 99.9%, APS<40μm
NRE-11019	CESIUM CARBONATE NANOPATICLES	Cs ₂ CO ₃ , 99.9%, APS<40μm
NRE-11021	CESIUM TUNGSTEN MICRO POWDER	Cs ₂ WO ₄ , 99.9%, APS<40μm
NRE-11022	CHROMIUM NITRIDE MICRO POWDER	CrN, 99.9%, APS<40μm
NRE-11023	COPPER MONOSULFIDE MICRO POWDER	CuS, 99.9%, APS<40μm
NRE-11024	GALLIUM SULFATE MICRO POWDER	Ga ₂ (SO ₄), 99.9%, APS<40μm
NRE-11025	GALLIUM TRICHLORIDE MICRO POWDER	GaCl ₃ , 99.9%, APS<40μm
NRE-11026	GOLD CHLORIDE MICRO POWDER	AuCl, 99.9%, APS<40μm
NRE-11027	HAFNIUM CARBIDE MICRO POWDER	HfC, 99.9%, APS<40μm
NRE-11028	HAFNIUM CHLORIDE MICRO POWDER	HfCl ₄ , 99.9%, APS<40μm
NRE-11029	HAFNIUM DIBORIDE MICRO POWDER	HfB ₂ , 99.9%, APS<40μm
NRE-11030	HAFNIUM SILICIDE MICRO POWDER	HfSi, 99.9%, APS<40μm
NRE-11031	INDIUM CHLORIDE MICRO POWDER	InCl ₃ , 99.9%, APS<40μm
NRE-11032	INDIUM FLUORIDE MICRO POWDER	InF ₃ , 99.9%, APS<40μm
NRE-11033	INDIUM NITRIDE MICRO POWDER	In ₃ N ₃ ·XH ₂ O, 99.9%, APS<40μm
NRE-11034	INDIUM SULFATE MICRO POWDER	In ₂ (SO ₄) ₃ , 99.9%, APS<40μm
NRE-11035	INDIUM SULFIDE MICRO POWDER	In ₂ S ₃ , 99.9%, APS<40μm
NRE-11036	INDIUM TELLURIDE NANOPATICLES	In ₂ Te ₃ , 99.9%, APS<40μm
NRE-11037	IRIDIUM CHLORIDE MICRO POWDER	IrCl ₃ , 99.9%, APS<40μm
NRE-11038	IRON BROMIDE MICRO POWDER	FeBr ₂ , 99.9%, APS<40μm
NRE-11039	IRON SULFIDE MICRO POWDER	FeS, 99.9%, APS<40μm
NRE-11040	LANTHANUM CHLORIDE MICRO POWDER	LaCl ₃ , 99.9%, APS<40μm
NRE-11041	LANTHANUM HEXABORIDE MICRO POWDER	LaB ₆ , 99.9%, APS<40μm
NRE-11042	LEAD NITRIDE MICRO POWDER	Pb(NO ₃) ₂ , 99.9%, APS<40μm
NRE-11043	LEAD TELLURIDE MICRO POWDER	PbTe, 99.9%, APS<40μm
NRE-11044	LEAD TUNGSTEN MICRO POWDER	PbWO ₄ , 99.9%, APS<40μm
NRE-11045	LITHIUM BOROHYDRIDE MICRO POWDER	LiBH ₄ , 99.9%, APS<40μm
NRE-11046	LITHIUM TELLURIDE MICRO POWDER	Li ₂ Te, 99.9%, APS<40μm

CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-11047	LITHIUM TUNGSTEN MICRO POWDER	Li ₂ WO ₄ , 99.9%, APS<40μm
NRE-11048	MAGNESIUM ALUMINATE MICRO POWDER	Al ₂ MgO ₄ , 99.9%, APS<40μm
NRE-11049	MAGNESIUM CHLORIDE MICRO POWDER	MgCl ₂ , 99.9%, APS<40μm
NRE-11050	MAGNESIUM FLURIDE MICRO POWDER	MgF ₂ , 99.9%, APS<40μm
NRE-11051	MAGNESIUM NITRIDE MICRO POWDER	Mg ₃ N ₂ , 99.9%, APS<40μm
NRE-11052	MAGNESIUM TUNGSTEN NANOPATICLES	MgWO ₄ , 99.9%, APS<40μm
NRE-11053	MANGANESE NITRIDE MICRO POWDER	Mn ₃ N ₂ , 99.9%, APS<40μm
NRE-11054	MANGANESE TUNGSTEN MICRO POWDER	MnO ₄ W, 99.9%, APS<40μm
NRE-11055	MOLYBDENUM CARBIDE MICRO POWDER	Mo ₂ C, 99.9%, APS<40μm
NRE-11056	MOLYBDENUM DISILICIDE MICRO POWDER	MoSi ₂ , 99.9%, APS<40μm
NRE-11057	MOLYBDENUM DISULFIDE MICRO POWDER	MoS ₂ , 99.9%, APS<40μm
NRE-11058	NIOBIUM CARBIDE MICRO POWDER	NbC, 99.9%, APS<40μm
NRE-11059	POTASSIUM TELLURITE MICRO POWDER	99.9%, APS<40μm
NRE-11060	POTASSIUM TITANATE MICRO POWDER	K ₂ O.8TiO ₂ , 99.9%, APS<40μm
NRE-11061	POTASSIUM TUNGSTATE MICRO POWDER	K ₂ WO ₄ , 99.9%, APS<40μm
NRE-11062	RUBIDIUM TUNGSTATE MICRO POWDER	Rb ₂ WO ₄ , 99.9%, APS<40μm
NRE-11063	SILICON CARBIDE MICRO POWDER	SiC, 99.9%, APS<40μm
NRE-11064	SILICON NITRIDE MICRO POWDER	Si ₃ N ₄ , 99.9%, APS<40μm
NRE-11065	SILVER TUNGSTATE MICRO POWDER	Ag ₂ WO ₄ , 99.9%, APS<40μm
NRE-11066	SODIUM FLURIDE MICRO POWDER	NaF, 99.9%, APS<40μm
NRE-11067	STRONTIUM FLUORIDE MICRO POWDER	SrF ₂ , 99.9%, APS<40μm
NRE-11068	STRONTIUM TUNGSTATE MICRO POWDER	SrWO ₄ , 99.9%, APS<40μm
NRE-11069	SUPER GRADE MICA MICRO POWDER	99.9%, APS<40μm
NRE-11070	TANTALUM NIOBIUM CARBIDE MICRO POWDER	TaCNbC, 99.9%, APS<40μm
NRE-11071	TELLURIUM CHLORIDE MICRO POWDER	Cl ₄ Te, 99.9%, APS<40μm
NRE-11072	THALLIUM CARBONATE MICRO POWDER	Tl ₂ CO ₃ , 99.9%, APS<40μm
NRE-11073	THALLIUM NITRATE MICRO POWDER	TlNO ₃ , 99.9%, APS<40μm
NRE-11074	THALLIUM SULFATE MICRO POWDER	Tl ₂ SO ₄ , 99.9%, APS<40μm
NRE-11075	TIN CHLORIDE MICRO POWDER	SnCl ₄ , 99.9%, APS<40μm
NRE-11076	TIN SELENIDE MICRO POWDER	SnSe, 99.9%, APS<40μm
NRE-11077	TITANIUM ALUMINIUM CARBIDE MICRO POWDER	Ti ₃ AlC ₂ , 99.9%, APS<40μm
NRE-11078	TITANIUM BORIDE MICRO POWDER	TiB ₂ , 99.9%, APS<40μm
NRE-11079	TITANIUM CARBIDE MICRO POWDER	TiC, 99.9%, APS<40μm
NRE-11080	TITANIUM DISILICIDE MICRO POWDER	TiSi ₂ , 99.9%, APS<40μm
NRE-11081	TITANIUM HYDRIDE MICRO POWDER	TiH ₂ , 99.9%, APS<40μm
NRE-11082	TITANIUM NITRIDE MICRO POWDER	TiN, 99.9%, APS<40μm
NRE-11083	TUNGSTEN CARBIDE COBALT MICRO POWDER	WC:CoC, 99.9%, APS<40μm
NRE-11084	TUNGSTEN CARBIDE MICRO POWDER	WC, 99.9%, APS<40μm
NRE-11085	TUNGSTEN TITANIUM CARBIDE MICRO POWDER	WC:TiC=50:50, 99.9%, APS<40μm
NRE-11086	TUNGSTEN TITANIUM TANTALUM CARBIDE MICRO POWDER	WC:TiC=50:50, 99.9%, APS<40μm
NRE-11087	VANADIUM SULFATE MICRO POWDER	V ₂ (SO ₄) ₃ , 99.9%, APS<40μm
NRE-11088	ZINC CHLORIDE MICRO POWDER	ZnCl ₂ , 99.9%, APS<40μm
NRE-11089	ZINC SULPHIDE MICRO POWDER	ZnS, 99.9%, APS<40μm
NRE-11090	ZINC TITANATE MICRO POWDER	ZnTiO ₃ , 99.9%, APS<40μm
NRE-11091	ZIRCONIUM CARBIDE MICRO POWDER	ZrC, 99.9%, APS<40μm

CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-11092	ZIRCONIUM DIBORIDE MICRO POWDER	ZrB ₂ , 99.9%, APS<40µm
NRE-11093	ZIRCONIUM DISILICIDE MICRO POWDER	ZrSi ₂ , 99.9%, APS<40µm
NRE-11094	ZIRCONIUM TUNGSTATE MICRO POWDER	Zr(WO ₄) ₂ , 99.9%, APS<40µm
NRE-11095	ZIRCONIUM(II) HYDRIDE MICRO POWDER	ZrH ₂ , 99.9%, APS<40µm
NRE-11096	ALUMINIUM NITRIDE POWDER	AlN, 99.9%, APS<40µm
NRE-11097	BARIUM CARBONATE POWDER	BaCO ₃ , 99.9%, APS<40µm
NRE-11098	BARIUM FERRITE POWDER	BaFe ₁₂ O ₁₉ , 99.9%, APS<40µm
NRE-11099	CALCIUM CARBONATE POWDER	CaCO ₃ , 99.9%, APS<40µm
NRE-11100	CESIUM TELLURIDE POWDER	Cs ₂ Te, 99.9%, APS<40µm
NRE-11101	COBALT SULPHIDE POWDER	CoS ₂ , 99.9%, APS<40µm
NRE-11102	COBALT ZINC FERRITE POWDER	CoFe ₂ O ₄ , 99.9%, APS<40µm
NRE-11103	COPPER SILICATE POWDER	CuO ₃ Si, 99.9%, APS<40µm
NRE-11104	HYDRATED LIME POWDER	Ca(OH) ₂ , 99.9%, APS<40µm
NRE-11105	IRON CARBIDE POWDER	Fe ₃ C, 99.9%, APS<40µm
NRE-11106	LANTHANUM FLUORIDE POWDER	LaF ₃ , 99.9%, APS<40µm
NRE-11107	LANTHANUM TRIFLUORIDE POWDER	LaF ₃ , 99.9%, APS<40µm
NRE-11108	LEAD SULFIDE POWDER	PbS, 99.9%, APS<40µm
NRE-11109	LITHIUM FLUORIDE POWDER	LiF, 99.9%, APS<40µm
NRE-11110	MAGNESIUM CARBONATE POWDER	MgCO ₃ , 99.9%, APS<40µm
NRE-11111	MAGNESIUM NITRATE POWDER	Mg(NO ₃) ₂ , 99.9%, APS<40µm
NRE-11112	MANGANESE SULFIDE POWDER	MnS, 99.9%, APS<40µm
NRE-11113	NICKEL ZINC FERRITE POWDER	NiZnFe ₂ O ₄ , 99.9%, APS<40µm
NRE-11114	POTASSIUM SULFATE POWDER	K ₂ SO ₄ , 99.9%, APS<40µm
NRE-11115	PTFE POWDER	99.9%, APS<40µm
NRE-11116	SILVER FLUORIDE POWDER	AgF, 99.9%, APS<40µm
NRE-11117	TANTALUM CARBIDE POWDER	TaC, 99.9%, APS<40µm
NRE-11118	TIN II SULFIDE POWDER	SnS, 99.9%, APS<40µm
NRE-11119	TIN NITRIDE POWDER	SnN, 99.9%, APS<40µm
NRE-11120	TITANIUM DIBORIDE POWDER	TiB ₂ , 99.9%, APS<40µm
NRE-11121	TUNGSTEN DISULFIDE POWDER	WS ₂ , 99.9%, APS<40µm
NRE-11122	VANADIUM CARBIDE POWDER	VC, 99.9%, APS<40µm
NRE-11123	ZINC CARBONATE POWDER	ZnCO ₃ , 99.9%, APS<40µm
NRE-11124	ZIRCONIUM DIHYDRIDE POWDER	ZrH ₂ , 99.9%, APS<40µm
NRE-11125	ZIRCONIUM NITRIDE POWDER	ZrN, 99.9%, APS<40µm
NRE-11126	ZIRCONIUM TUNGSTATE POWDER	Zr(WO ₄) ₂ , 99.9%, APS<40µm
NRE-11127	VANADIUM HYDRIDE POWDER	H ₅ V-5, 99.9%, APS<40µm
NRE-11128	HAFNIUM HYDRIDE POWDER	HfH ₂ , 99.9%, APS<40µm
NRE-11129	CALCIUM HYDRIDE POWDER	CaH ₂ , 99.9%, APS<40µm
NRE-11130	ALUMINIUM CARBIDE POWDER	Al ₄ C ₃ , 99.9%, APS<40µm
NRE-11131	MANGANESE CARBIDE POWDER	C ₂ Mg, 99.9%, APS<40µm
NRE-11132	TITANIUM CARBONITRIDE POWDER	TiCN, 99.9%, APS<40µm
NRE-11133	CHROMIUM BORIDE POWDER	Cr ₃ B ₄ , 99.9%, APS<40µm
NRE-11134	NICKEL BORIDE POWDER	Ni ₂ B, 99.9%, APS<40µm
NRE-11135	MAGNESIUM BORIDE POWDER	MgB ₂ , 99.9%, APS<40µm
NRE-11136	SILICON BORIDE POWDER	SiB ₆ , 99.9%, APS<40µm

CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-11137	ALUMINIUM BORIDE POWDER	AlB ₂ , 99.9%, APS<40μm
NRE-11138	MOLYBDENUM BORIDE POWDER	BMo, 99.9%, APS<40μm
NRE-11139	MANGANESE BORIDE POWDER	MnB, 99.9%, APS<40μm
NRE-11140	NIOBIUM BORIDE POWDER	NbB ₂ , 99.9%, APS<40μm
NRE-11141	COBALT BORIDE POWDER	Co ₂ B, 99.9%, APS<40μm
NRE-11142	VANADIUM BORIDE POWDER	B ₂ V, 99.9%, APS<40μm
NRE-11143	TUNGSTEN BORIDE POWDER	B ₅ W ₂ , 99.9%, APS<40μm
NRE-11144	CALCIUM NITRIDE POWDER	Ca ₃ N ₂ , 99.9%, APS<40μm
NRE-11145	HAFNIUM NITRIDE POWDER	HfN, 99.9%, APS<40μm
NRE-11146	ALUMINIUM SILICIDE POWDER	Al ₄ Si ₃ , 99.9%, APS<40μm
NRE-11147	CHROMIUM SILICIDE POWDER	CrSi ₂ , 99.9%, APS<40μm
NRE-11148	MAGNESIUM SILICIDE POWDER	Mg ₂ Si, 99.9%, APS<40μm
NRE-11149	CALCIUM SILICIDE POWDER	CaSi ₂ , 99.9%, APS<40μm
NRE-11150	VANADIUM SILICIDE POWDER	VS ₂ , 99.9%, APS<40μm
NRE-11151	MANGANESE SILICIDE POWDER	MnSi ₂ , 99.9%, APS<40μm
NRE-11152	FERRIC SILICIDE POWDER	FeSi, 99.9%, APS<40μm
NRE-11153	NIOBIUM SILICIDE POWDER	NbSi ₂ , 99.9%, APS<40μm
NRE-11154	TUNGSTEN SILICIDE POWDER	WS ₂ , 99.9%, APS<40μm
NRE-11155	NICKEL SILICIDE POWDER	Ni ₂ Si, 99.9%, APS<40μm




CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-57001	CELLULOSE POWDER	99%+, APS <40 µm
NRE-57002	CHITOSAN POWDER	99%+, APS <40 µm
NRE-57003	CLAY POWDER	99%+, APS <40 µm
NRE-57004	EXPANDED PERLITE POWDER	99%+, APS <40 µm
NRE-57005	GELATIN POWDER	99%+, APS <40 µm
NRE-57006	HALLOYSITE CLAY	99%+, APS <40 µm
NRE-57007	HYDROXYAPATITE POWDER	99%+, APS <40 µm
NRE-57008	KAOLIN CLAY	99%+, APS <40 µm
NRE-57009	MICA POWDER	99%+, APS <40 µm
NRE-57010	MONTMORILLONITE BENTONITE CLAY	99%+, APS <40 µm
NRE-57011	MONTMORILLONITE CLAY	99%+, APS <40 µm
NRE-57012	CLAY MODIFIED ASPHALT MATERIALS	99%+, APS <40 µm
NRE-57013	CLAY POWDER	99%+, APS <40 µm
NRE-57014	PERLITE CLAY MIX POWDER	99%+, APS <40 µm
NRE-57015	PERLITE CLAY POWDER	99%+, APS <40 µm
NRE-57016	SAPONITE CLAY	99%+, APS <40 µm
NRE-57017	SILVER COATED HYDROXYAPATITE POWDER	99%+, APS <40 µm
NRE-57018	ZEOLITE POWDER	99%+, APS <40 µm
NRE-57019	SUPER GRADE MICA POWDER	99%+, APS <40 µm
NRE-57020	CALCINED MICA POWDER	99%+, APS <40 µm
NRE-57021	BENTONITE CLAY	99%+, APS <40 µm




CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-12001	ALUMINUM NANORODS	PURITY 99%+
NRE-12002	ANTIMONY NANORODS	PURITY 99%+
NRE-12003	ARSENIC NANORODS	PURITY 99%+
NRE-12004	BISMUTH NANORODS	PURITY 99%+
NRE-12005	BORON NANORODS	PURITY 99%+
NRE-12006	CADMIUM NANORODS	PURITY 99%+
NRE-12007	CARBON NANORODS	PURITY 99%+
NRE-12008	CHROMIUM NANORODS	PURITY 99%+
NRE-12009	COBALT NANORODS	PURITY 99%+
NRE-12010	COPPER NANORODS	PURITY 99%+
NRE-12011	GERMANIUM NANORODS	PURITY 99%+
NRE-12012	GOLD NANORODS	PURITY 99%+
NRE-12013	GOLD NANORODS, PALLADIUM-COATED	PURITY 99%+
NRE-12014	GOLD NANORODS, PLATINUM-COATED	PURITY 99%+
NRE-12015	HAFNIUM NANORODS	PURITY 99%+
NRE-12016	INDIUM NANORODS	PURITY 99%+
NRE-12017	IRIDIUM NANORODS	PURITY 99%+
NRE-12018	IRON NANORODS	PURITY 99%+
NRE-12019	LEAD NANORODS	PURITY 99%+
NRE-12020	MANGANESE NANORODS	PURITY 99%+
NRE-12021	MANGANESE(IV) OXIDE NANORODS	PURITY 99%+
NRE-12022	MOLYBDENUM NANORODS	PURITY 99%+
NRE-12023	NICKEL NANORODS	PURITY 99%+
NRE-12024	NIOBIUM NANORODS	PURITY 99%+
NRE-12025	OSMIUM NANORODS	PURITY 99%+
NRE-12026	PALLADIUM NANORODS	PURITY 99%+
NRE-12027	PLATINUM NANORODS	PURITY 99%+
NRE-12028	RHENIUM NANORODS	PURITY 99%+
NRE-12029	RHODIUM NANORODS	PURITY 99%+
NRE-12030	RUTHENIUM NANORODS	PURITY 99%+
NRE-12031	SCANDIUM NANORODS	PURITY 99%+
NRE-12032	SILICON NANORODS	PURITY 99%+
NRE-12033	SILVER NANORODS	PURITY 99%+
NRE-12034	TANTALUM NANORODS	PURITY 99%+
NRE-12035	THALLIUM NANORODS	PURITY 99%+
NRE-12036	TIN NANORODS	PURITY 99%+
NRE-12037	TITANIUM NANORODS	PURITY 99%+
NRE-12038	TUNGSTEN NANORODS	PURITY 99%+
NRE-12039	VANADIUM NANORODS	PURITY 99%+
NRE-12040	YTTRIUM NANORODS	PURITY 99%+
NRE-12041	ZINC NANORODS	PURITY 99%+
NRE-12042	ZIRCONIUM NANORODS	PURITY 99%+



CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-13001	ALUMINUM OXIDE NANOWIRES	PURITY 99%+
NRE-13002	COPPER NANOWIRES	PURITY 99%+
NRE-13003	GERMANIUM NANOWIRES	PURITY 99%+
NRE-13004	GOLD NANOWIRES	PURITY 99%+
NRE-13005	LEAD NANOWIRES	PURITY 99%+
NRE-13006	NICKEL NANOWIRES	PURITY 99%+
NRE-13007	NICKEL OXIDE NANOFIBERS	PURITY 99%+
NRE-13008	NICKEL(II) OXIDE NANOWIRES	PURITY 99%+
NRE-13009	SILICON NANOWIRES	PURITY 99%+
NRE-13010	SILVER NANOWIRES	PURITY 99%+
NRE-13011	SODIUM AMMONIUM TRIMOLYBDATE NANOWIRES	PURITY 99%+
NRE-13012	TIN OXIDE NANOWIRES	PURITY 99%+
NRE-13013	TITANIUM(IV) OXIDE NANOWIRE	PURITY 99%+
NRE-13014	TUNGSTEN(VI) OXIDE NANOWIRES	PURITY 99%+
NRE-13015	ZINC OXIDE NANOWIRE	PURITY 99%+



CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-14001	ALUMINUM CARBON NANOTUBES	PURITY 99%+
NRE-14002	ALUMINUM NITRIDE NANOTUBES	PURITY 99%+
NRE-14003	BORON NITRIDE NANOTUBES	PURITY 99%+
NRE-14004	CADMIUM SELENIDE NANOTUBE	PURITY 99%+
NRE-14005	COPPER CARBON NANOTUBES	PURITY 99%+
NRE-14006	COPPER NANOTUBES	PURITY 99%+
NRE-14007	DOUBLE-WALLED CARBON NANOTUBES	PURITY 99%+
NRE-14008	FULLERENE NANOTUBES, MULTI-WALLED	PURITY 99%+
NRE-14009	FULLERENE NANOTUBES, SINGLE-WALLED	PURITY 99%+
NRE-14010	GOLD NANOTUBES	PURITY 99%+
NRE-14011	HALLOYSITE NANOTUBES	PURITY 99%+
NRE-14012	MAGNESIUM CARBON NANOTUBES	PURITY 99%+
NRE-14013	MOLYBDENUM DISULFIDE NANOTUBES	PURITY 99%+
NRE-14014	MULTI-WALLED CARBON NANOTUBES	PURITY 99%+
NRE-14015	NICKEL CARBON NANOTUBES	PURITY 99%+
NRE-14016	SILICON NANOTUBES	PURITY 99%+
NRE-14017	SILVER GOLD NANOTUBE	PURITY 99%+
NRE-14018	SILVER NANOTUBES	PURITY 99%+
NRE-14019	TIN CARBON NANOTUBES	PURITY 99%+
NRE-14020	TITANIUM CARBON NANOTUBES	PURITY 99%+
NRE-14021	TITANIUM OXIDE NANOTUBES	PURITY 99%+





CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-15001	CADMIUM SELENIDE HOLLOW NANOSPHERES	99%+, APS <100 nm
NRE-15002	COPPER(I) OXIDE NANOSPHERES	99%+, APS <100 nm
NRE-15003	LEAD SULFIDE HOLLOW NANOSPHERES	99%+, APS <100 nm
NRE-15004	PLATINUM HOLLOW NANOSPHERES	99%+, APS <100 nm
NRE-15005	SILICA MICROSPHERES	99%+, APS <100 nm
NRE-15006	SILICON DIOXIDE NANOSPHERES	99%+, APS <100 nm
NRE-15007	SILICON OXIDE HOLLOW NANOSPHERES	99%+, APS <100 nm
NRE-15008	SILVER GOLD HOLLOW NANOSPHERES	99%+, APS <100 nm
NRE-15009	TITANIUM DIOXIDE NANOSPHERES	99%+, APS <100 nm




CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-16001	ALUMINIUM SILICON OXIDE CORE SHELL NANOPARTICLES	Al/SiO ₂ , 99.9%, APS<100nm
NRE-16002	CADMIUM SELENIDE CADMIUM SULFUR SILICA CORE SHELL NP'S	CdSe/CdS/SiO ₂ , 99.9%, APS<100nm
NRE-16003	CADMIUM SELENIDE SILICA CORE SHELL NANOPARTICLES	CdSe/SiO ₂ , 99.9%, APS<100nm
NRE-16004	CADMIUM SELENIUM ZINC SULFUR SILICA CORE SHELL NP'S	CdSe/ ZnS/SiO ₂ , 99.9%, APS<100nm
NRE-16005	CADMIUM SELENIUM, CADMIUM SULFIDE/ SILICON CORE SHELL NP'S	CdSe/CdS/SiO ₂ , 99.9%, APS<100nm
NRE-16006	CADMIUM SELENIUM/ SILICA CORE SHELL NANOPARTICLES	CdSe/SiO ₂ , 99.9%, APS<100nm
NRE-16007	CADMIUM SELENIUM/CADMIUM SULFIDE CORE SHELL NP'S	CdSe CdS, 99.9%, APS<100nm
NRE-16008	CADMIUM SULFIDE TITANIUM OXIDE CORE SHELL NP'S	CdS/TiO ₂ , 99.9%, APS<100nm
NRE-16009	CADMIUM SULFIDE/ MERCURY SULFIDE CORE SHELL NP'S	CdS/HgS, 99.9%, APS<100nm
NRE-16010	CADMIUM SULFIDE/ SILVER CORE SHELL NANOPARTICLES	CdS/Ag, 99.9%, APS<100nm
NRE-16011	CADMIUM SULFIDE/CADMIUM SELENIUM CORE SHELL NP'S	CdS/CdSe, 99.9%, APS<100nm
NRE-16012	CADMIUM SULFIDE/MERCURY SULFIDE/ CADMIUM SULFIDE CORE SHELL	CdS/HgS, 99.9%, APS<100nm
NRE-16013	CADMIUM SULFIDE/ZINC SULFIDE CORE SHELL NP'S	CdS/ZnS, 99.9%, APS<100nm
NRE-16014	CADMIUM SULFUR SILICA CORE SHELL NANOPARTICLES	CdS/SiO ₂ , 99.9%, APS<100nm
NRE-16015	CADMIUM TELLURIDE/ SILICON CORE SHELL NANOPARTICLES	CdTe/SiO ₂ , 99.9%, APS<100nm
NRE-16016	CALCIUM COBALT SILICA CORE SHELL NANOPARTICLES	CaCo ₃ /SiO ₂ , 99.9%, APS<100nm
NRE-16017	CALCIUM OXIDE/IRON OXIDE CORE SHELL NANOPARTICLES	CaO/Fe ₂ O ₃ , 99.9%, APS<100nm
NRE-16018	CARBON NANOTUBES/ TIN OXIDE CORE SHELL NANOPARTICLES	CNT/SnO ₂ , 99.9%, APS<100nm
NRE-16019	CARBON/ GOLD CORE SHELL NANOPARTICLES	C/Au, 99.9%, APS<100nm
NRE-16020	CdS Ag ₂ S CORE SHELL NANOPARTICLES	CdS Ag ₂ S, 99.9%, APS<100nm
NRE-16021	CdS/ PbS CORE SHELL NANOPARTICLES	CdS/PbS, 99.9%, APS<100nm
NRE-16022	CADMIUM SULPHIDE/SILVER SULFIDE CORE SHELL NP'S	CdS/Ag ₂ S, 99.9%, APS<100nm
NRE-16023	CADMIUM SELENIDE/CADMIUM SULPHIDE/ZINC SULPHIDE CORE SHELL	CdSe CdS ZnS, 99.9%, APS<100nm
NRE-16024	CdSe ZnSe CORE SHELL NANOPARTICLES	CdSe ZnSe, 99.9%, APS<100nm
NRE-16025	CdSe ZnTe CORE SHELL NANOPARTICLES	CdSe ZnTe, 99.9%, APS<100nm
NRE-16026	CdSe/ZnS CORE SHELL NANOPARTICLES	CdSe/ZnS, 99.9%, APS<100nm
NRE-16027	CADMIUM TELLURIUM/CADMIUM SULFIDE CORE SHELL NP'S	CdTe CdS, 99.9%, APS<100nm
NRE-16028	COBALT CADMIUM SELENIUM CORE SHELL NANOPARTICLES	Co/CdSe, 99.9%, APS<100nm
NRE-16029	COBALT SILICA CORE SHELL NANOPARTICLES	Co/SiO ₂ , 99.9%, APS<100nm
NRE-16030	COPPER OXIDE/ ALUMINUM CORE SHELL NANOPARTICLES	Cu ₂ O/Al, 99.9%, APS<100nm
NRE-16031	COPPER/ COPPER OXIDE CORE SHELL NANOPARTICLES	Cu/Cu ₂ O, 99.9%, APS<100nm
NRE-16032	ZINC SULFIDE/ SILICON OXIDE CORE SHELL NANO PARTICLES	99.9%, APS<100nm
NRE-16033	IRON OXIDE SILICA CORE SHELL NANOPARTICLES	Fe ₃ O ₄ SiO ₂ , 99.9%, APS<100nm
NRE-16034	GADOLINIUM OXIDE THULIUM SILICA CORE SHELL NP'S	99.9%, APS<100nm
NRE-16035	GOLD CADMIUM SELENIDE CORE SHELL NANOPARTICLES	Au/CdSe, 99.9%, APS<100nm
NRE-16036	GOLD CADMIUM SELENIUM CORE SHELL NANOPARTICLES	Au/CdSe, 99.9%, APS<100nm
NRE-16037	GOLD CARBON CORE SHELL NANOPARTICLES	Au/C, 99.9%, APS<100nm
NRE-16038	GOLD COBALT CORE SHELL NANOPARTICLES	Au/Co, 99.9%, APS<100nm
NRE-16039	GOLD COPPER IODIDE CORE SHELL NANOPARTICLES	Au/CuI, 99.9%, APS<100nm
NRE-16040	GOLD IRON OXIDE CORE SHELL NANOPARTICLES	Au/Fe ₂ O ₃ , 99.9%, APS<100nm
NRE-16041	GOLD PALLADIUM CORE SHELL NANOPARTICLES	Au/Pd, 99.9%, APS<100nm
NRE-16042	GOLD PLATINUM CORE SHELL NANOPARTICLES	Au/Pt, 99.9%, APS<100nm
NRE-16043	GOLD SILICA CORE SHELL NANOPARTICLES	Au/SiO ₂ , 99.9%, APS<100nm
NRE-16044	GOLD SILICON OXIDE CORE-SHELL NANOPARTICLES	Au/SiO ₂ , 99.9%, APS<100nm
NRE-16045	GOLD SILVER CORE SHELL NANOPARTICLES	Au/ Ag, 99.9%, APS<100nm

CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-16046	GOLD TITANIUM DIOXIDE CORE SHELL NANOPARTICLES	Au/TiO ₂ , 99.9%, APS<100nm
NRE-16047	IRON CARBON CORE SHELL NANOPARTICLES	Fe/C, 99.9%, APS<100nm
NRE-16048	IRON COPPER GOLD PLATINUM PALLADIUM SILVER CORE SHELL NP'S	Fe, Cu/Au, Pd, Pt, Cu, 99.9%, APS<100nm
NRE-16049	IRON NICKEL SILICA CORE SHELL NANOPARTICLES	Fe Ni/SiO ₂ , 99.9%, APS<100nm
NRE-16050	IRON OXIDE DIETHYLAMINE ETHYL CORE SHELL NANOPARTICLES	Fe ₂ O ₃ /C ₄ H ₁₁ N, 99.9%, APS<100nm
NRE-16051	IRON OXIDE GOLD CORE SHELL NANOPARTICLES	Fe ₂ O ₃ /Au, 99.9%, APS<100nm
NRE-16052	IRON OXIDE POLYVINYL CORE SHELL NANOPARTICLES	Fe ₂ O ₃ /C ₄ H ₁₁ N, 99.9%, APS<100nm
NRE-16053	IRON OXIDE POLYVINYL ALCOHOL CORE SHELL NANOPARTICLES	Fe ₂ O ₃ / (C ₂ H ₄ O) _x , 99.9%, APS<100nm
NRE-16054	IRON OXIDE SILICA CORE SHELL NANOPARTICLES	Fe/SiO ₂ , 99.9%, APS<100nm
NRE-16055	IRON OXIDE/ CARBON CORE SHELL NANOPARTICLES	Fe ₃ O ₄ /C, 99.9%, APS<100nm
NRE-16056	IRON OXIDE/ GOLD CORE SHELL NANO PARTICLES	Fe ₃ O ₄ /Au, 99.9%, APS<100nm
NRE-16057	IRON OXIDE/ SILICON OXIDE/ ALUMINIUM OXIDE CORE SHELL NP'S	Fe ₃ O ₄ /SiO ₂ Al ₂ O ₃ , 99.9%, APS<100nm
NRE-16058	IRON OXIDE/TITANIUM OXIDE CORE SHELL NANOPARTICLES	Fe ₃ O ₄ /TiO ₂ , 99.9%, APS<100nm
NRE-16059	IRON PLATINUM CADMIUM SULFUR CORE SHELL NANOPARTICLES	Fe Pt/CdS, 99.9%, APS<100nm
NRE-16060	IRON PLATINUM CORE SHELL NANOPARTICLES	Fe/Pt, 99.9%, APS<100nm
NRE-16061	IRON PLATINUM/ IRON OXIDE CORE SHELL NANOPARTICLES	Fe ₅₈ Pt ₄₂ /Fe ₃ O ₄ , 99.9%, APS<100nm
NRE-16062	IRON SILVER CORE SHELL NANOPARTICLES	Fe/Ag, 99.9%, APS<100nm
NRE-16063	LANTHANUM FLUORINE EUROPIUM LANTHANUM FLUORINE	LaF ₃ / Eu _{0.2} La _{0.8} F ₃ , 99.9%, APS<100nm
NRE-16064	LANTHANUM POLONIUM EUROPIUM/ LANTHANUM POLONIUM	LaPo ₄ : Eu ³⁺ / LaPo ₄ , 99.9%, APS<100nm
NRE-16065	LANTHANUM STRONTIUM MANGANESE GOLD CORE SHELL	La ₂ /3Sr1/3MnO ₃ /Au, 99.9%, APS<100nm
NRE-16066	MAGNESIUM OXIDE CORE SHELL NANOPARTICLES	MgO Fe ₂ O ₃ , 99.9%, APS<100nm
NRE-16067	NICKEL CARBON CORE SHELL NANOPARTICLES	Ni/C, 99.9%, APS<100nm
NRE-16068	NICKEL PLATINUM CORE SHELL NANOPARTICLES	Ni/Pt, 99.9%, APS<100nm
NRE-16069	NICKEL SILICA CORE SHELL NANOPARTICLES	Ni/SiO ₂ , 99.9%, APS<100nm
NRE-16070	NICKEL SILVER CORE SHELL NANOPARTICLES	Ni/Ag, 99.9%, APS<100nm
NRE-16071	PALLADIUM PALLADIUM OXIDE CORE SHELL NANOPARTICLES	Pd/PdO, 99.9%, APS<100nm
NRE-16072	LEAD SELENIUM/CADMIUM SELENIUM CORE SHELL NANOPARTICLES	PbTe/CdSe, 99.9%, APS<100nm
NRE-16073	ZINC CADMIUM SELENIUM/ZINC SULFIDE CORE SHELL NANOPARTICLES	ZnCdSe/ZnS, 99.9%, APS<100nm
NRE-16074	SILICA LATEX RUBBER CORE-SHELL NANOPARTICLES	SiO ₂ , 99.9%, APS<100nm
NRE-16075	SILICON OXIDE/ GOLD CORE SHELL NANOPARTICLES	SiO ₂ /Au, 99.9%, APS<100nm
NRE-16076	SILICON OXIDE/ ZINC SILICON OXIDE MANGANESE CORE SHELL NP'S	SiO ₂ /Zn ₂ SiO ₄ :Mn ₂ , 99.9%, APS<100nm
NRE-16077	SILVER CARBON CORE SHELL NANOPARTICLES	Ag/C, 99.9%, APS<100nm
NRE-16078	SILVER INDIUM SILICON OXIDE CORE SHELL NANOPARTICLES	AgIn/SiO ₂ , 99.9%, APS<100nm
NRE-16079	SILVER IODIDE SILICON OXIDE CORE SHELL NANOPARTICLES	AgI/SiO ₂ , 99.9%, APS<100nm
NRE-16080	SILVER SILICA CORE SHELL NANOPARTICLES	Ag/SiO ₂ , 99.9%, APS<100nm
NRE-16081	SILVER SILICON OXIDE CORE SHELL NANOPARTICLE	Ag/SiO ₂ , 99.9%, APS<100nm
NRE-16082	SILVER SILVER SELENIDE CORE SHELL NANOPARTICLES	Ag/Ag ₂ Se, 99.9%, APS<100nm
NRE-16083	YTTERBIUM SILICA CORE SHELL NANOPARTICLES	Yb ³⁺ /SiO ₂ , 99.9%, APS<100nm
NRE-16084	ZINC OXIDE/ SILVER CORE SHELL NANOPARTICLES	ZnO/Ag, 99.9%, APS<100nm
NRE-16085	ZINC OXIDE/ TITANIUM OXIDE CORE SHELL NANOPARTICLES	ZnO/TiO ₂ , 99.9%, APS<100nm
NRE-16086	ZINC SILICON OXIDE EUROPIUM SILICON OXIDE CORE SHELL NP'S	Zn ₂ SiO ₄ :Eu ³⁺ /SiO ₂ , 99.9%, APS<100nm
NRE-16087	ZINC SULFIDE MANGANESE/ ZINC OXIDE CORE SHELL NP'S	ZnS:Mn/ZnO, 99.9%, APS<100nm
NRE-16088	ZINC SULFIDE/CADMIUM SELENIUM CORE SHELL NANO PARTICLES	ZnS/CdSe, 99.9%, APS<100nm
NRE-16089	ZINC SULFUR MANGANESE SILICA CORE SHELL NANOPARTICLES	ZnS:Mn/SiO ₂ , 99.9%, APS<100nm
NRE-16090	ZINC/ ZINC OXIDE CORE SHELL NANOPARTICLES	Zn/ ZnO, 99.9%, APS<100nm

CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-16091	ZINC OXIDE/SILICA CORE SHELL NANOPARTICLES	ZnO SiO ₂ , 99.9%, APS<100nm
NRE-16092	ZINC SELENIUM/CADMIUM SELENIUM CORE SHELL NANOPARTICLES	ZnSe CdSe, 99.9%, APS<100nm
NRE-16093	COPPER SILVER COATED NANOPARTICLES	CuAg, 99.9%, APS<100nm
NRE-16094	CADMIUM TELLURIUM SILICA CORE SHELL NANOPARTICLES	CdTe/SiO ₂ , 99.9%, APS<100nm
NRE-16095	CADMIUM TELLURIUM CADMIUM SELENIUM CORE SHELL NP'S	CdTe/CdSe, 99.9%, APS<100nm



CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-17001	ALUMINIUM METAL ORGANIC FRAMEWORK	Al-MIL-53, 99%
NRE-17002	CHROMIUM METAL ORGANIC FRAMEWORK	Cr-MIL-100, 99%
NRE-17003	CHROMIUM MOF	Cr-MIL-101, 99%
NRE-17004	COPPER METAL ORGANIC FRAMEWORK	Cu (HKUST-1), 99%
NRE-17005	DYSPROSIUM METAL ORGANIC FRAMEWORK	Dy-MOF, 99%
NRE-17006	EUROPIUM METAL ORGANIC FRAMEWORK	Eu-MOF-76, 99%
NRE-17007	FERROUS METAL ORGANIC FRAMEWORK	Fe-MOF, 99%
NRE-17008	IRON METAL ORGANIC FRAMEWORK	Fe-(MOF-5), 99%
NRE-17009	MANGANESE METAL ORGANIC FRAMEWORK	Mn, 99%
NRE-17010	NANOPOROUS IRON METAL ORGANIC FRAMEWORKS	MIL-89, 99%
NRE-17011	NICKEL METAL ORGANIC FRAMEWORK	Ni-MOF, 99%
NRE-17012	SAMARIUM METAL ORGANIC FRAMEWORK	Sm-MOFS, 99%
NRE-17013	TERBIUM METAL ORGANIC FRAMEWORK	Tb-MOF-76, 99%
NRE-17014	TIN METAL ORGANIC FRAMEWORK	Sn-MOF, 99%
NRE-17015	TITANIUM METAL ORGANIC FRAMEWORK	Ti-MOF, 99%
NRE-17016	VANADIUM METAL ORGANIC FRAMEWORKS	V-MIL-47, 99%
NRE-17017	ZINC METAL ORGANIC FRAMEWORK	Zn-ZIF-8, 99%
NRE-17018	ZIRCONIUM METAL ORGANIC FRAMEWORK	Zr-MOF, 99%



CODE NO.	PRODUCT NAME
NRE-18001	CdSeS/ZnS ALLOYED QD COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 490 nm, 6 nm DIAMETER, 1 mg/mL IN H ₂ O
NRE-18002	CdSeS/ZnS ALLOYED QD COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 525 nm, 6 nm DIAMETER, 1 mg/mL IN H ₂ O
NRE-18003	CdSeS/ZnS ALLOYED QD COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 575 nm, 6 nm DIAMETER, 1 mg/mL IN H ₂ O
NRE-18004	CdSeS/ZnS ALLOYED QD COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 630 nm, 6 nm DIAMETER, 1 mg/mL IN H ₂ O
NRE-18005	CdSeS/ZnS ALLOYED QD COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 665 nm, 6 nm DIAMETER, 1 mg/mL IN H ₂ O
NRE-18006	CdSeS/ZnS ALLOYED QD FLUORESCENCE λ_{em} 450 nm, 6 nm DIAMETER, 1 mg/mL IN TOLUENE
NRE-18007	CdSeS/ZnS ALLOYED QD FLUORESCENCE λ_{em} 490 nm, 6 nm DIAMETER, 1 mg/mL IN TOLUENE
NRE-18008	CdSeS/ZnS ALLOYED QD FLUORESCENCE λ_{em} 525 nm, 6 nm DIAMETER, 1 mg/mL IN TOLUENE
NRE-18009	CdSeS/ZnS ALLOYED QD FLUORESCENCE λ_{em} 540 nm, 6 nm DIAMETER, 1 mg/mL IN H ₂ O
NRE-18010	CdSeS/ZnS ALLOYED QD FLUORESCENCE λ_{em} 540 nm, 6 nm DIAMETER, 1 mg/mL IN TOLUENE
NRE-18011	CdSeS/ZnS ALLOYED QD FLUORESCENCE λ_{em} 575 nm, 6 nm DIAMETER, 1 mg/mL IN TOLUENE
NRE-18012	CdSeS/ZnS ALLOYED QD FLUORESCENCE λ_{em} 630 nm, 6 nm DIAMETER, 1 mg/mL IN TOLUENE
NRE-18013	CdSeS/ZnS ALLOYED QD FLUORESCENCE λ_{em} 665 nm, 6 nm DIAMETER, 1 mg/mL IN TOLUENE





CODE NO.	PRODUCT NAME
NRE-19001	InP/ZnS QD STABILIZED WITH OLEYLAMINE LIGANDS, FLUORESCENCE λ_{em} 530 nm, 5 mg/mL IN TOLUENE
NRE-19002	InP/ZnS QD STABILIZED WITH OLEYLAMINE LIGANDS, FLUORESCENCE λ_{em} 560 nm, 5 mg/mL IN TOLUENE
NRE-19003	InP/ZnS QD STABILIZED WITH OLEYLAMINE LIGANDS, FLUORESCENCE λ_{em} 590 nm, 5 mg/mL IN TOLUENE
NRE-19004	InP/ZnS QD STABILIZED WITH OLEYLAMINE LIGANDS, FLUORESCENCE λ_{em} 620 nm, 5 mg/mL IN TOLUENE
NRE-19005	InP/ZnS QD STABILIZED WITH OLEYLAMINE LIGANDS, FLUORESCENCE λ_{em} 650 nm, 5 mg/mL IN TOLUENE
NRE-19006	PbS CORE-TYPE QD OLEIC ACID COATED, FLUORESCENCE λ_{em} 1000 nm, 10 mg/mL IN TOLUENE
NRE-19007	PbS CORE-TYPE QD OLEIC ACID COATED, FLUORESCENCE λ_{em} 1200 nm, 10 mg/mL IN TOLUENE
NRE-19008	PbS CORE-TYPE QD OLEIC ACID COATED, FLUORESCENCE λ_{em} 1400 nm, 10 mg/mL IN TOLUENE
NRE-19009	PbS CORE-TYPE QD OLEIC ACID COATED, FLUORESCENCE λ_{em} 1600 nm, 10 mg/mL IN TOLUENE



CODE NO.	PRODUCT NAME
NRE-20001	CdTe CORE-TYPE QUANTUM DOTS COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 510 nm, POWDER
NRE-20002	CdTe CORE-TYPE QUANTUM DOTS COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 520 nm, POWDER
NRE-20003	CdTe CORE-TYPE QUANTUM DOTS COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 570 nm, POWDER
NRE-20004	CdTe CORE-TYPE QUANTUM DOTS COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 610 nm, POWDER
NRE-20005	CdTe CORE-TYPE QUANTUM DOTS COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 710 nm, POWDER
NRE-20006	CdTe CORE-TYPE QUANTUM DOTS COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 770 nm, POWDER
NRE-20007	PbS CORE-TYPE QUANTUM DOTS OLEIC ACID COATED, FLUORESCENCE λ_{em} 1000 nm, 10 mg/mL IN TOLUENE
NRE-20008	PbS CORE-TYPE QUANTUM DOTS OLEIC ACID COATED, FLUORESCENCE λ_{em} 1200 nm, 10 mg/mL IN TOLUENE
NRE-20009	PbS CORE-TYPE QUANTUM DOTS OLEIC ACID COATED, FLUORESCENCE λ_{em} 1400 nm, 10 mg/mL IN TOLUENE
NRE-20010	PbS CORE-TYPE QUANTUM DOTS OLEIC ACID COATED, FLUORESCENCE λ_{em} 1600 nm, 10 mg/mL IN TOLUENE
NRE-20011	CdS CORE-TYPE QUANTUM DOTS, 5 mg/mL IN TOLUENE
NRE-20012	CdSeS/ZnS ALLOYED QD COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 490 nm
NRE-20013	CdSeS/ZnS ALLOYED QUANTUM DOTS COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 525 nm
NRE-20014	CdSeS/ZnS ALLOYED QUANTUM DOTS COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 575 nm
NRE-20015	CdSeS/ZnS ALLOYED QUANTUM DOTS COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 630 nm
NRE-20016	CdSeS/ZnS ALLOYED QUANTUM DOTS COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 665 nm
NRE-20017	CdSeS/ZnS ALLOYED QUANTUM DOTS FLUORESCENCE λ_{em} 450 nm, 6 nm DIAMETER
NRE-20018	CdSeS/ZnS ALLOYED QUANTUM DOTS FLUORESCENCE λ_{em} 490 nm, 6 nm DIAMETER
NRE-20019	CdSeS/ZnS ALLOYED QUANTUM DOTS FLUORESCENCE λ_{em} 525 nm, 6 nm DIAMETER
NRE-20020	CdSeS/ZnS ALLOYED QUANTUM DOTS FLUORESCENCE λ_{em} 540 nm, 6 nm DIAMETER
NRE-20021	CdSeS/ZnS ALLOYED QUANTUM DOTS FLUORESCENCE λ_{em} 540 nm, 6 nm DIAMETER
NRE-20022	CdSeS/ZnS ALLOYED QUANTUM DOTS FLUORESCENCE λ_{em} 575 nm, 6 nm DIAMETER
NRE-20023	CdSeS/ZnS ALLOYED QUANTUM DOTS FLUORESCENCE λ_{em} 630 nm, 6 nm DIAMETER
NRE-20024	CdSeS/ZnS ALLOYED QUANTUM DOTS KIT FLUORESCENCE λ_{em} 490-665 nm, 6 nm DIAMETER


Quantum Dots



CODE NO.	PRODUCT NAME
NRE-21001	CdSe/ZnS CORE SHELL (CdSe/ZINC SULPHIDE QD-520nm)
NRE-21002	CdS/ZnS CORE-SHELL TYPE QD LYOPHILIZED, λ_{em} 400 nm, SOLID
NRE-21003	CdS/ZnS CORE-SHELL TYPE QD LYOPHILIZED, λ_{em} 450 nm, SOLID
NRE-21004	CdS/ZnS CORE-SHELL TYPE QD OLEIC ACID FUNCTIONALIZED, λ_{em} 425 nm, 5 mg/mL IN TOLUENE
NRE-21005	CdS/ZnS CORE-SHELL TYPE QD OLEIC ACID FUNCTIONALIZED, λ_{em} 450 nm, 5 mg/mL IN TOLUENE
NRE-21006	CdS/ZnS CORE-SHELL TYPE QD STABILIZED WITH OLEYLAMINE LIGANDS, λ_{em} 530 nm, 5 mg/mL IN TOLUENE
NRE-21007	CdS/ZnS CORE-SHELL TYPE QD STABILIZED WITH OLEYLAMINE LIGANDS, λ_{em} 590 nm, 5 mg/mL IN TOLUENE
NRE-21008	CdSe ZnS CORE SHELL QUANTUM DOT (CADMIUM SELENIDE/ZINC SULPHIDE QD-500nm)
NRE-21009	CdSe/CdS CORE-SHELL TYPE QUANTUM RODS λ_{em} 530 nm, 5 mg/mL IN HEXANE
NRE-21010	CdSe/CdS CORE-SHELL TYPE QUANTUM RODS λ_{em} 560 nm, 5 mg/mL IN HEXANE
NRE-21011	CdSe/CdS CORE-SHELL TYPE QUANTUM RODS λ_{em} 590 nm, 5 mg/mL IN HEXANE
NRE-21012	CdSe/CdS CORE-SHELL TYPE QUANTUM RODS λ_{em} 620 nm, 5 mg/mL IN HEXANE
NRE-21013	CdSe/ZnS 510, CORE-SHELL TYPE QD, 5 mg/mL IN TOLUENE
NRE-21014	CdSe/ZnS 560, CORE-SHELL TYPE QD, 5 mg/mL IN TOLUENE
NRE-21015	CdSe/ZnS 590, CORE-SHELL TYPE QD, 5 mg/mL IN TOLUENE
NRE-21016	CdSe/ZnS CORE-SHELL TYPE QD AMINE FUNCTIONALIZED, λ_{em} 400 nm, 4 μ m in 10 mm PBS
NRE-21017	CdSe/ZnS CORE-SHELL TYPE QD AMINE FUNCTIONALIZED, λ_{em} 425 nm, 4 μ m in 10 mm PBS
NRE-21018	CdSe/ZnS CORE-SHELL TYPE QD AMINE FUNCTIONALIZED, λ_{em} 450 nm, 4 μ m in 10 mm PBS
NRE-21019	CdSe/ZnS CORE-SHELL TYPE QD AMINE FUNCTIONALIZED, λ_{em} 520 nm, 4 μ m in 10 mm PBS
NRE-21020	CdSe/ZnS CORE-SHELL TYPE QD AMINE FUNCTIONALIZED, λ_{em} 540 nm, 4 μ m in 10 mm PBS
NRE-21021	CdSe/ZnS CORE-SHELL TYPE QD AMINE FUNCTIONALIZED, λ_{em} 560 nm, 4 μ m in 10 mm PBS
NRE-21022	CdSe/ZnS CORE-SHELL TYPE QD AMINE FUNCTIONALIZED, λ_{em} 580 nm, 4 μ m in 10 mm PBS
NRE-21023	CdSe/ZnS CORE-SHELL TYPE QD AMINE FUNCTIONALIZED, λ_{em} 620 nm, 4 μ m in 10 mm PBS
NRE-21024	CdSe/ZnS CORE-SHELL TYPE QD AMINE FUNCTIONALIZED, λ_{em} 630 nm, 4 μ m in 10 mm PBS
NRE-21025	CdSe/ZnS CORE-SHELL TYPE QD AMINE FUNCTIONALIZED, λ_{em} 645 nm, 4 μ m in 10 mm PBS
NRE-21026	CdSe/ZnS CORE-SHELL TYPE QD AMINE FUNCTIONALIZED, λ_{em} 665 nm, 4 μ m in 10 mm PBS
NRE-21027	CdSe/ZnS CORE-SHELL TYPE QD AMINE FUNCTIONALIZED, λ_{em} 600 nm, 4 μ m in 10 mm PBS
NRE-21028	CdSe/ZnS CORE-SHELL TYPE QD CARBOXYLIC ACID FUNCTIONALIZED, λ_{em} 400 nm, 1 mg/mL IN H2O
NRE-21029	CdSe/ZnS CORE-SHELL TYPE QD CARBOXYLIC ACID FUNCTIONALIZED, λ_{em} 425 nm, 1 mg/mL IN H2O
NRE-21030	CdSe/ZnS CORE-SHELL TYPE QD CARBOXYLIC ACID FUNCTIONALIZED, λ_{em} 450 nm, 1 mg/mL IN H2O
NRE-21031	CdSe/ZnS CORE-SHELL TYPE QD CARBOXYLIC ACID FUNCTIONALIZED, λ_{em} 540 nm, 1 mg/mL IN H2O
NRE-21032	CdSe/ZnS CORE-SHELL TYPE QD CARBOXYLIC ACID FUNCTIONALIZED, λ_{em} 560 nm, 1 mg/mL IN H2O
NRE-21033	CdSe/ZnS CORE-SHELL TYPE QD CARBOXYLIC ACID FUNCTIONALIZED, λ_{em} 580 nm, 1 mg/mL IN H2O
NRE-21034	CdSe/ZnS CORE-SHELL TYPE QD CARBOXYLIC ACID FUNCTIONALIZED, λ_{em} 620 nm, 1 mg/mL IN H2O
NRE-21035	CdSe/ZnS CORE-SHELL TYPE QD CARBOXYLIC ACID FUNCTIONALIZED, λ_{em} 630 nm, 1 mg/mL IN H2O
NRE-21036	CdSe/ZnS CORE-SHELL TYPE QD CARBOXYLIC ACID FUNCTIONALIZED, λ_{em} 645 nm, 1 mg/mL IN H2O
NRE-21037	CdSe/ZnS CORE-SHELL TYPE QD CARBOXYLIC ACID FUNCTIONALIZED, λ_{em} 665 nm, 1 mg/mL IN H2O
NRE-21038	CdSe/ZnS CORE-SHELL TYPE QD CARBOXYLIC ACID FUNCTIONALIZED, λ_{em} 600 nm, 1 mg/mL IN H2O
NRE-21039	CdSe/ZnS CORE-SHELL TYPE QD LYOPHILIZED, λ_{em} 425 nm, SOLID
NRE-21040	CdSe/ZnS CORE-SHELL TYPE QD LYOPHILIZED, λ_{em} 540 nm, SOLID
NRE-21041	CdSe/ZnS CORE-SHELL TYPE QD LYOPHILIZED, λ_{em} 560 nm, SOLID
NRE-21042	CdSe/ZnS CORE-SHELL TYPE QD LYOPHILIZED, λ_{em} 580 nm, SOLID
NRE-21043	CdSe/ZnS CORE-SHELL TYPE QD LYOPHILIZED, λ_{em} 600 nm, SOLID
NRE-21044	CdSe/ZnS CORE-SHELL TYPE QD LYOPHILIZED, λ_{em} 620 nm, SOLID
NRE-21045	CdSe/ZnS CORE-SHELL TYPE QD LYOPHILIZED, λ_{em} 630 nm, SOLID

CODE NO.	PRODUCT NAME
NRE-21046	CdSe/ZnS CORE-SHELL TYPE QD LYOPHILIZED, λ_{em} 665 nm, SOLID
NRE-21047	CdSe/ZnS CORE-SHELL TYPE QD OLEIC ACID FUNCTIONALIZED, λ_{em} 400 nm, 5 mg/mL IN TOLUENE
NRE-21048	CdSe/ZnS CORE-SHELL TYPE QD OLEIC ACID FUNCTIONALIZED, λ_{em} 425 nm, SOLID
NRE-21049	CdSe/ZnS CORE-SHELL TYPE QD OLEIC ACID FUNCTIONALIZED, λ_{em} 450 nm, SOLID
NRE-21050	CdSe/ZnS CORE-SHELL TYPE QD PDDA COATED, λ_{em} 400 nm, 4 μ m IN H ₂ O
NRE-21051	CdSe/ZnS CORE-SHELL TYPE QD PDDA COATED, λ_{em} 425 nm, 4 μ m IN H ₂ O
NRE-21052	CdSe/ZnS CORE-SHELL TYPE QD PDDA COATED, λ_{em} 450 nm, 4 μ m IN H ₂ O
NRE-21053	CdSe/ZnS CORE-SHELL TYPE QD PDDA COATED, λ_{em} 520 nm, 4 μ m IN H ₂ O
NRE-21054	CdSe/ZnS CORE-SHELL TYPE QD PDDA COATED, λ_{em} 540 nm, 4 μ m IN H ₂ O
NRE-21055	CdSe/ZnS CORE-SHELL TYPE QD PDDA COATED, λ_{em} 560 nm, 4 μ m IN H ₂ O
NRE-21056	CdSe/ZnS CORE-SHELL TYPE QD PDDA COATED, λ_{em} 580 nm, 4 μ m IN H ₂ O
NRE-21057	CdSe/ZnS CORE-SHELL TYPE QD PDDA COATED, λ_{em} 620nm, 4 μ m IN H ₂ O
NRE-21058	CdSe/ZnS CORE-SHELL TYPE QD PDDA COATED, λ_{em} 630nm, 4 μ m IN H ₂ O
NRE-21059	CdSe/ZnS CORE-SHELL TYPE QD PDDA COATED, λ_{em} 645 nm, 4 μ m IN H ₂ O
NRE-21060	CdSe/ZnS CORE-SHELL TYPE QD PEG FUNCTIONALIZED, λ_{em} 400 nm, 4 μ m IN H ₂ O
NRE-21061	CdSe/ZnS CORE-SHELL TYPE QD PEG FUNCTIONALIZED, λ_{em} 425 nm, 4 μ m IN H ₂ O
NRE-21062	CdSe/ZnS CORE-SHELL TYPE QD PEG FUNCTIONALIZED, λ_{em} 450 nm, 4 μ m IN H ₂ O
NRE-21063	CdSe/ZnS CORE-SHELL TYPE QD PEG FUNCTIONALIZED, λ_{em} 520 nm, 4 μ m IN H ₂ O
NRE-21064	CdSe/ZnS CORE-SHELL TYPE QD PEG FUNCTIONALIZED, λ_{em} 540 nm, 4 μ m IN H ₂ O
NRE-21065	CdSe/ZnS CORE-SHELL TYPE QD PEG FUNCTIONALIZED, λ_{em} 560 nm, 4 μ m IN H ₂ O
NRE-21066	CdSe/ZnS CORE-SHELL TYPE QD PEG FUNCTIONALIZED, λ_{em} 580 nm, 4 μ m IN H ₂ O
NRE-21067	CdSe/ZnS CORE-SHELL TYPE QD PEG FUNCTIONALIZED, λ_{em} 600 nm, 4 μ m IN H ₂ O
NRE-21068	CdSe/ZnS CORE-SHELL TYPE QD PEG FUNCTIONALIZED, λ_{em} 620 nm, 4 μ m IN H ₂ O
NRE-21069	CdSe/ZnS CORE-SHELL TYPE QD PEG FUNCTIONALIZED, λ_{em} 645 nm, 4 μ m IN H ₂ O
NRE-21070	CdSe/ZnS CORE-SHELL TYPE QD STABILIZED WITH OCTADECYLAMINE LIGANDS, λ_{em} 520 nm, 5 mg/mL IN TOLUENE
NRE-21071	CdSe/ZnS CORE-SHELL TYPE QD STABILIZED WITH OCTADECYLAMINE LIGANDS, λ_{em} 520 nm, SOLID
NRE-21072	CdSe/ZnS CORE-SHELL TYPE QD STABILIZED WITH OCTADECYLAMINE LIGANDS, λ_{em} 540 nm, 5 mg/mL IN TOLUENE
NRE-21073	CdSe/ZnS CORE-SHELL TYPE QD STABILIZED WITH OCTADECYLAMINE LIGANDS, λ_{em} 540nm, SOLID
NRE-21074	CdSe/ZnS CORE-SHELL TYPE QD STABILIZED WITH OCTADECYLAMINE LIGANDS, λ_{em} 560 nm, 5 mg/mL IN TOLUENE
NRE-21075	CdSe/ZnS CORE-SHELL TYPE QD STABILIZED WITH OCTADECYLAMINE LIGANDS, λ_{em} 560 nm, SOLID
NRE-21076	CdSe/ZnS CORE-SHELL TYPE QD STABILIZED WITH OCTADECYLAMINE LIGANDS, λ_{em} 580 nm, 5 mg/mL IN TOLUENE
NRE-21077	CdSe/ZnS CORE-SHELL TYPE QD STABILIZED WITH OCTADECYLAMINE LIGANDS, λ_{em} 580 nm, SOLID
NRE-21078	CdSe/ZnS CORE-SHELL TYPE QD STABILIZED WITH OCTADECYLAMINE LIGANDS, λ_{em} 600 nm, 5 mg/mL IN TOLUENE
NRE-21079	CdSe/ZnS CORE-SHELL TYPE QD STABILIZED WITH OCTADECYLAMINE LIGANDS, λ_{em} 600 nm, SOLID
NRE-21080	CdSe/ZnS CORE-SHELL TYPE QD STABILIZED WITH OCTADECYLAMINE LIGANDS, λ_{em} 620 nm, SOLID
NRE-21081	CdSe/ZnS CORE-SHELL TYPE QD STABILIZED WITH OCTADECYLAMINE LIGANDS, λ_{em} 630 nm, SOLID
NRE-21082	CdSe/ZnS CORE-SHELL TYPE QD STABILIZED WITH OCTADECYLAMINE LIGANDS, λ_{em} 645 nm, 5 mg/mL IN TOLUENE
NRE-21083	CdSe/ZnS CORE-SHELL TYPE QD STABILIZED WITH OCTADECYLAMINE LIGANDS, λ_{em} 645 nm, SOLID
NRE-21084	CdSe/ZnS CORE-SHELL TYPE QUANTUM DOTS STABILIZED WITH OCTADECYLAMINE LIGANDS, FLUORESCENCE λ_{em} 665nm
NRE-21085	CdSeS/ZnS ALLOYED QUANTUM DOTS COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 490 nm
NRE-21086	CdSeS/ZnS ALLOYED QUANTUM DOTS COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 525nm
NRE-21087	CdSeS/ZnS ALLOYED QUANTUM DOTS COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 540 nm
NRE-21088	CdSeS/ZnS ALLOYED QUANTUM DOTS COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 575 nm
NRE-21089	CdSeS/ZnS ALLOYED QUANTUM DOTS COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 630 nm
NRE-21090	CdSeS/ZnS ALLOYED QUANTUM DOTS COOH FUNCTIONALIZED, FLUORESCENCE λ_{em} 665nm

CODE NO.	PRODUCT NAME
NRE-21091	CdSeS/ZnS ALLOYED QUANTUM DOTS FLUORESCENCE λ_{em} 450 nm, 6 nm DIAMETER
NRE-21092	CdSeS/ZnS ALLOYED QUANTUM DOTS FLUORESCENCE λ_{em} 490 nm, 6 nm DIAMETER
NRE-21093	CdSeS/ZnS ALLOYED QUANTUM DOTS FLUORESCENCE λ_{em} 525 nm, 6 nm DIAMETER
NRE-21094	CdSeS/ZnS ALLOYED QUANTUM DOTS FLUORESCENCE λ_{em} 540 nm, 6 nm DIAMETER
NRE-21095	CdSeS/ZnS ALLOYED QUANTUM DOTS FLUORESCENCE λ_{em} 575 nm, 6 nm DIAMETER
NRE-21096	CdSeS/ZnS ALLOYED QUANTUM DOTS FLUORESCENCE λ_{em} 630 nm, 6 nm DIAMETER
NRE-21097	CdSeS/ZnS ALLOYED QUANTUM DOTS FLUORESCENCE λ_{em} 665 nm, 6 nm DIAMETER
NRE-21098	CdSeS/ZnS ALLOYED QUANTUM DOTS KIT FLUORESCENCE λ_{em} 490-665 nm, 6 nm DIAMETER
NRE-21099	InP/ZnS QUANTUM DOTS STABILIZED WITH OLEYLAMINE LIGANDS, FLUORESCENCE λ_{em} 560 nm
NRE-21100	InP/ZnS QUANTUM DOTS STABILIZED WITH OLEYLAMINE LIGANDS, FLUORESCENCE λ_{em} 620 nm
NRE-21101	InP/ZnS QUANTUM DOTS STABILIZED WITH OLEYLAMINE LIGANDS, FLUORESCENCE λ_{em} 650 nm



CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-22001	COPPER NANODISPERSION SOLVENT-WATER	Cu, 99.9 %, APS<50nm, 50PPM
NRE-22002	COPPER NANODISPERSION SOLVENT-WATER	Cu, 99.9 %, APS<50nm, 100PPM
NRE-22003	COPPER NANODISPERSION SOLVENT-WATER	Cu, 99.9 %, APS 50nm, 500PPM
NRE-22004	COPPER NANODISPERSION SOLVENT-WATER	Cu, 99.9 %, APS 50nm, 1000PPM
NRE-22005	GOLD NANODIPERSION SOLVENT-WATER	Au, 99.9%, APS<20 nm, 50PPM
NRE-22006	GOLD NANODIPERSION SOLVENT-WATER	Au, 99.9%, APS<20 nm, 100PPM
NRE-22007	GOLD NANODIPERSION SOLVENT-WATER	Au, 99.9%, APS<20 nm, 500PPM
NRE-22008	GOLD NANODIPERSION SOLVENT-WATER	Au, 99.9%, APS 20 nm, 1000PPM
NRE-22009	GOLD NANODIPERSION SOLVENT-WATER	Au, 99.9%, APS<20 nm, 50PPM
NRE-22010	GOLD NANODIPERSION SOLVENT-ETHANOL	Au, 99.9%, APS<20 nm, 100PPM
NRE-22011	GOLD NANODIPERSION SOLVENT-ETHANOL	Au, 99.9%, APS<20 nm, 500PPM
NRE-22012	GOLD NANODIPERSION SOLVENT-ETHANOL	Au, 99.9%, APS<20 nm, 1000PPM
NRE-22013	MAGNESIUM NANODISPERSION	Mg, 99.9 %, APS<80nm
NRE-22014	NIOBIUM NANODISPERSION WITH EDA LIGAND	99.9%, APS 80-100nm
NRE-22015	PLATINUM NANODIPERSION SOLVENT-WATER	Pt, 99.9 %, APS<50nm, 50PPM
NRE-22016	PLATINUM NANODIPERSION SOLVENT-WATER	Pt, 99.9 %, APS<50nm, 100PPM
NRE-22017	PLATINUM NANODIPERSION SOLVENT-WATER	Pt, 99.9 %, APS<50nm, 500PPM
NRE-22018	PLATINUM NANODIPERSION SOLVENT-WATER	Pt, 99.9 %, APS<50nm, 1000PPM
NRE-22019	PLATINUM NANODIPERSION SOLVENT-ETHANOL	Pt, 99.9 %, APS<50nm, 50PPM
NRE-22020	PLATINUM NANODIPERSION SOLVENT-ETHANOL	Pt, 99.9 %, APS<50nm, 100PPM
NRE-22021	PLATINUM NANODIPERSION SOLVENT-ETHANOL	Pt, 99.9 %, APS<50nm, 500PPM
NRE-22022	PLATINUM NANODIPERSION SOLVENT-ETHANOL	Pt, 99.9 %, APS<50nm, 1000PPM
NRE-22023	SILVER NANODIPERSION SOLVENT-WATER	Ag, 99.9 %, APS<20nm, 50PPM
NRE-22024	SILVER NANODIPERSION SOLVENT-WATER	Ag, 99.9 %, APS<20nm, 100PPM
NRE-22025	SILVER NANODIPERSION SOLVENT-WATER	Ag, 99.9 %, APS<20nm, 500PPM
NRE-22026	SILVER NANODIPERSION SOLVENT-WATER	Ag, 99.9 %, APS<20nm, 1000PPM
NRE-22027	SILVER NANODIPERSION SOLVENT-ETHANOL	Ag, 99.9 %, APS<20nm, 50PPM
NRE-22028	SILVER NANODIPERSION SOLVENT-ETHANOL	Ag, 99.9 %, APS<20nm, 100PPM
NRE-22029	SILVER NANODIPERSION SOLVENT-ETHANOL	Ag, 99.9 %, APS<20nm, 500PPM
NRE-22030	SILVER NANODIPERSION SOLVENT-ETHANOL	Ag, 99.9 %, APS<20nm, 1000PPM
NRE-22031	NANO SILVER FLUORIDE CHITOSAN STABILIZED	Ag, 99.9% APS <20nm, 200PPM
NRE-22032	SILVER NANODISPERSION SOLVENT-WATER	Ag, 99.9% APS <2nm, 200PPM
NRE-22033	SILVER NANODISPERSION SOLVENT-WATER	Ag, 99.9% APS <15nm, 5000PPM
NRE-22034	SILVER NANODISPERSION SOLVENT-WATER	Ag, 99.9% APS <15nm, 50000PPM
NRE-22035	SILVER NANODISPERSION SOLVENT-WATER	Ag, 99.9% APS <2nm, 2000PPM
NRE-22036	SILVER CONDUCTIVE PASTE	Ag <60%, APS <2μm)
NRE-22037	SILVER NANODIPERSION SOLVENT-ETHANOL	Ag, 99.9 %, APS <15nm, 500-5000PPM
NRE-22038	SILVER NANODIPERSION SOLVENT-ETHANOL	Ag, 99.9 %, APS <2nm, 100-10000PPM
NRE-22039	DIAMOND NANODISPERSION SOLVENT-ETHANOL	C, 98.3 %, APS <10nm
NRE-22040	DIAMOND NANODISPERSION WATER-ETHANOL	C, 98.3 %, APS <10nm


CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-23001	ALUMINUM HYDROXIDE NANODISPERSION SOLVENT-WATER	Al(OH) ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23002	ALUMINUM HYDROXIDE NANODISPERSION SOLVENT-ETHANOL	Al(OH) ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23003	ALUMINUM OXIDE NANODISPERSION SOLVENT-ETHYLENE GLYCOL	Al ₂ O ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23004	ALUMINUM OXIDE NANODISPERSION SOLVENT-WATER	Al ₂ O ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23005	ANTIMONY TIN OXIDE NANODISPERSION SOLVENT-WATER	SnO ₂ :Sb ₂ O ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23006	ANTIMONY TIN OXIDE NANODISPERSION SOLVENT-ETHANOL	SnO ₂ :Sb ₂ O ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23007	ALUMINIUM DOPED ZINC OXIDE (AZO) NANODISPERSION SOL-WATER	Al-ZnO, 99.9%, APS<50 nm, 1000PPM
NRE-23008	ALUMINIUM DOPED ZINC OXIDE (AZO) NANODISPERSION SOL-ETHANOL	Al-ZnO, 99.9%, APS<50 nm, 1000PPM
NRE-23009	BARIUM TITANATE NANODISPERSION SOLVENT-WATER	BaTiO ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23010	BARIUM TITANATE NANODISPERSION SOLVENT-ETHANOL	BaTiO ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23011	BISMUTH OXIDE NANODISPERSION SOLVENT-WATER	Bi ₂ O ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23012	BISMUTH OXIDE NANODISPERSION SOLVENT-ETHANOL	Bi ₂ O ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23013	CERIUM OXIDE NANODISPERSION SOLVENT-WATER	CeO ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23014	CERIUM OXIDE NANODISPERSION SOLVENT-ETHANOL	CeO ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23015	CHROMIUM OXIDE NANODISPERSION SOLVENT-WATER	Cr ₂ O ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23016	CHROMIUM OXIDE NANODISPERSION SOLVENT-ETHANOL	Cr ₂ O ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23017	COBALT OXIDE NANODISPERSION SOLVENT-WATER	Co ₃ O ₄ , 99.9%, APS<50 nm, 1000PPM
NRE-23018	COBALT OXIDE NANODISPERSION SOLVENT-ETHANOL	Co ₃ O ₄ , 99.9%, APS<50 nm, 1000PPM
NRE-23019	COBALT OXIDE NANODISPERSION SOLVENT-ETHYLENE GLYCOL	Co ₃ O ₄ , 99.9%, APS<50 nm, 1000PPM
NRE-23020	COPPER OXIDE NANODISPERSION SOLVENT-WATER	CuO, 99.9%, APS<50 nm, 1000PPM
NRE-23021	COPPER OXIDE NANODISPERSION SOLVENT-ETHANOL	CuO, 99.9%, APS<50 nm, 1000PPM
NRE-23022	Fe(III) OXIDE HYDROXIDE NANODISPERSION SOLVENT-WATER	FeOOH, 99.9%, APS<50 nm, 1000PPM
NRE-23023	Fe(III) OXIDE HYDROXIDE NANODISPERSION SOLVENT-ETHANOL	FeOOH, 99.9%, APS<50 nm, 1000PPM
NRE-23024	Fe(III) OXIDE HYDROXIDE NANODISPERSION SOL-ETHYLENE GLYCOL	FeOOH, 99.9%, APS<50 nm, 1000PPM
NRE-23025	IRON OXIDE NANODISPERSION SOLVENT-WATER	Fe ₂ O ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23026	IRON OXIDE NANODISPERSION SOLVENT-ETHANOL	Fe ₂ O ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23027	IRON OXIDE NANODISPERSION SOLVENT-ETHYLENE GLYCOL	Fe ₂ O ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23028	IRON OXIDE NANODISPERSION SOLVENT-WATER	Fe ₃ O ₄ , 99.9%, APS<50 nm, 1000PPM
NRE-23029	IRON OXIDE NANODISPERSION SOLVENT-ETHANOL	Fe ₃ O ₄ , 99.9%, APS<50 nm, 1000PPM
NRE-23030	IRON OXIDE NANODISPERSION SOLVENT-ETHYLENE GLYCOL	Fe ₃ O ₄ , 99.9%, APS<50 nm, 1000PPM
NRE-23031	IRON HYDROXIDE NANODISPERSION SOLVENT-WATER	Fe(OH) ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23032	IRON HYDROXIDE NANODISPERSION SOLVENT-ETHANOL	Fe(OH) ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23033	IRON HYDROXIDE NANODISPERSION SOLVENT-ETHYLENE GLYCOL	Fe(OH) ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23034	INDIUM TIN OXIDE (ITO) NANODISPERSION SOLVENT-WATER	In ₂ O ₃ :SnO ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23035	INDIUM TIN OXIDE (ITO) NANODISPERSION SOLVENT-ETHANOL	In ₂ O ₃ :SnO ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23036	INDIUM TIN OXIDE (ITO) NANODISPERSION SOL-ETHYLENE GLYCOL	In ₂ O ₃ :SnO ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23037	MAGNESIUM HYDROXIDE NANODISPERSION SOLVENT-WATER	Mg(OH) ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23038	MAGNESIUM HYDROXIDE NANODISPERSION SOLVENT-ETHANOL	Mg(OH) ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23039	MAGNESIUM HYDROXIDE NANODISPERSION SOL-ETHYLENE GLYCOL	Mg(OH) ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23040	NICKEL HYDROXIDE NANODISPERSION SOLVENT-WATER	Ni(OH) ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23041	NICKEL HYDROXIDE NANODISPERSION SOLVENT-ETHANOL	Ni(OH) ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23042	NICKEL HYDROXIDE NANODISPERSION SOLVENT-ETHYLENE GLYCOL	Ni(OH) ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23043	CALCIUM CARBONATE NANODISPERSION SOLVENT-WATER	CaCO ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23044	CALCIUM CARBONATE NANODISPERSION SOLVENT-ETHANOL	CaCO ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23045	CALCIUM CARBONATE NANODISPERSION SOLVENT-ETHYLENE GLYCOL	CaCO ₃ , 99.9%, APS<50 nm, 1000PPM

CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-23046	SILICON OXIDE NANODISPERSION SOLVENT-WATER	SiO ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23047	SILICON OXIDE NANODISPERSION SOLVENT-ETHANOL	SiO ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23048	SILICON OXIDE NANODISPERSION SOLVENT-ETHYLENE GLYCOL	SiO ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23049	TIN DIOXIDE NANODISPERSION SOLVENT-WATER	SnO ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23050	TIN DIOXIDE NANODISPERSION SOLVENT-ETHANOL	SnO ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23051	TIN DIOXIDE NANODISPERSION SOLVENT-ETHYLENE GLYCOL	SnO ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23052	TITANIUM OXIDE NANODISPERSION SOLVENT-WATER	TiO ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23053	TITANIUM OXIDE NANODISPERSION SOLVENT-ETHANOL	TiO ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23054	TITANIUM OXIDE NANODISPERSION SOLVENT-ETHYLENE GLYCOL	TiO ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23055	TUNGSTEN OXIDE NANODISPERSION SOLVENT-WATER	WO ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23056	TUNGSTEN OXIDE NANODISPERSION SOLVENT-ETHANOL	WO ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23057	TUNGSTEN OXIDE NANODISPERSION SOLVENT-ETHYLENE GLYCOL	WO ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23058	YTTRIUM OXIDE NANODISPERSION SOLVENT-WATER	Y ₂ O ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23059	YTTRIUM OXIDE NANODISPERSION SOLVENT-ETHANOL	Y ₂ O ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23060	YTTRIUM OXIDE NANODISPERSION SOLVENT-ETHYLENE GLYCOL	Y ₂ O ₃ , 99.9%, APS<50 nm, 1000PPM
NRE-23061	ZINC OXIDE NANODISPERSION SOLVENT-WATER	ZnO, 99.9%, APS<50 nm, 1000PPM
NRE-23062	ZINC OXIDE NANODISPERSION SOLVENT-ETHANOL	ZnO, 99.9%, APS<50 nm, 1000PPM
NRE-23063	ZINC OXIDE NANODISPERSION SOLVENT-ETHYLENE GLYCOL	ZnO, 99.9%, APS<50 nm, 1000PPM
NRE-23064	ZIRCONIUM OXIDE NANODISPERSION SOLVENT-WATER	ZrO ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23065	ZIRCONIUM OXIDE NANODISPERSION SOLVENT-ETHANOL	ZrO ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23066	ZIRCONIUM OXIDE NANODISPERSION SOLVENT-ETHYLENE GLYCOL	ZrO ₂ , 99.9%, APS<50 nm, 1000PPM
NRE-23067	IRON HYDROXIDE NANODISPERSION SOLVENT-WATER	Fe(OH) ₃ , 99.9%, APS<5 nm





CODE NO.	PRODUCT NAME
NRE-24001	CONDUCTIVE CARBON BLACK NANOPARTICLES WATER DISPERSION
NRE-24002	FULLERENE FULLERENOLS 1000ppm DISPERSION
NRE-24003	GRAPHENE BASED CONDUCTIVE INKS (C, ORGANIC SOLVENT, FLAKE: 1-3nm)
NRE-24004	GRAPHENE CARBON NANOTUBES AG NANOPOWDER MIXED 6WT% IN WATER DISPERSION
NRE-24005	GRAPHENE CARBON NANOTUBES MIXED 6WT% IN WATER DISPERSION
NRE-24006	GRAPHENE DISPERSION ORGANIC SOLVENT (C, PURITY: >95 %, HYDROPHOBIC)
NRE-24007	GRAPHENE NANOPATELETS WATER DISPERSION (95+%, 6WT%, AVERAGE 3-6 LAYERS)
NRE-24008	HIGHLY CONDUCTIVE CARBON BLACK & GRAPHENE & CARBON NANOTUBES MIXED 20WT% IN WATER DISPERSION
NRE-24009	MWCNTs DISPERSION (>95%, DIAMETER: 20-30nm, LENGTH: 5-15µm)
NRE-24010	MWCNTs DOPED WITH NANOPOWDER Ag IN WATER DISPERSION, 15WT%
NRE-24011	MWCNTs ETHANOL DISPERSION (3WT%, >95%, OD: 20-30 nm, LENGTH 10-30µm)
NRE-24012	MWCNTs ETHANOL DISPERSION (3WT%, >95%, OD: 50-80 nm, LENGTH 10-20µm)
NRE-24013	MWCNTs ETHANOL DISPERSION (3WT%, >95+%, OD: 5-15 nm, LENGTH 50µm)
NRE-24014	MWCNTs ISOPROPANOL DISPERSION (3WT%, >95%, OD: 20-30 nm, LENGTH 10-30µm)
NRE-24015	MWCNTs ISOPROPANOL DISPERSION (3WT%, >95%, OD: 50-80 nm, LENGTH 10-20µm)
NRE-24016	MWCNTs ISOPROPANOL DISPERSION (3WT%, >95+%, OD: 5-15 nm, LENGTH 50µm)
NRE-24017	MWCNTs N-BUTANOL DISPERSION (3WT%, >95%, OD: 20-30 nm, LENGTH 10-30µm)
NRE-24018	MWCNTs N-BUTANOL DISPERSION (3WT%, >95%, OD: 50-80 nm, LENGTH 10-20µm)
NRE-24019	MWCNTs N-BUTANOL DISPERSION (3WT%, >95+%, OD: 5-15 NM, LENGTH 50µm)
NRE-24020	MWCNTs N-METHYL-2-PYRROLIDINONE DISPERSION (3WT%, >95%, OD: 20-30 nm, LENGTH 10-30µm)
NRE-24021	MWCNTs N-METHYL-2-PYRROLIDINONE DISPERSION (3WT%, >95%, OD: 50-80 nm, LENGTH 10-20µm)
NRE-24022	MWCNTs N-METHYL-2-PYRROLIDINONE DISPERSION (3WT%, >95+%, OD: 5-15 nm, LENGTH 50µm)
NRE-24023	MWCNTs WATER DISPERSION (15WT%, CNTS PURITY >95%, OD: 20-30 nm, LENGTH: 10-30µm)
NRE-24024	MWCNTs WATER DISPERSION (3WT%, >95%, OD: 20-30 nm, LENGTH 10-30µm)
NRE-24025	MWCNTs WATER DISPERSION (3WT%, >95%, OD: 50-80 nm, LENGTH 10-20µm)
NRE-24026	MWCNTs WATER DISPERSION (3WT%, >95+%, OD: 5-15 NM, LENGTH 50µm)
NRE-24027	POLYHYDROXYLATED FULLERENE (FULLERENOLS) / C60, -OH FUNCTIONALIZED
NRE-24028	PVP COATED GRAPHITIZED MWNTS / DRY SOLID DISPERSIBLE MWCNTs (>99.9%, OD: 10-20nm)
NRE-24029	PVP COATED MWCNTs / DRY SOLID DISPERSIBLE MWCNTs (>95%, OD: <7nm)
NRE-24030	PVP COATED MWCNTs / DRY SOLID DISPERSIBLE MWCNTs (>95%, OD: 20-30 nm)
NRE-24031	PVP COATED MWCNTs / DRY SOLID DISPERSIBLE MWCNTs (>95%, OD: 50-80 nm)
NRE-24032	PVP COATED SHORT SWCNTs / DRY SOLID DISPERSIBLE SHORT SWCNTs, PURITY >90%
NRE-24033	PVP COATED SWCNTs / DRY SOLID DISPERSIBLE SWCNTs, PURITY >90%
NRE-24034	RESEARCH GRADE GRAPHENE NANOPATELETS WATER DISPERSION
NRE-24035	RESEARCH GRADE SINGLE LAYER GRAPHENE NANOPARTICLES 1WT% WATER DISPERSION
NRE-24036	RESEARCH GRADE SINGLE LAYER GRAPHENE OXIDE WATER DISPERSION
NRE-24037	SHORT SWCNTs WATER DISPERSION
NRE-24038	SHORT SWCNTs DOPED WITH NANOPOWDER AG IN WATER DISPERSION, 1WT%
NRE-24039	SHORT SWCNTs ETHANOL DISPERSION
NRE-24040	SHORT SWCNTs XYLENE DISPERSION
NRE-24041	SHORT SWCNTs-COOH ETHANOL DISPERSION
NRE-24042	SHORT SWCNTs-COOH XYLENE DISPERSION
NRE-24043	SHORT SWCNTs-OCOH WATER DISPERSION
NRE-24044	SHORT SWCNTs-OH ETHANOL DISPERSION
NRE-24045	SHORT SWCNTs-OH WATERDISPERSION

CODE NO.	PRODUCT NAME
NRE-24046	SHORT SWCNTs-OH XYLENE DISPERSION
NRE-24047	SILICON GRAPHENE CARBON NANOTUBES MIXED 6WT% IN WATER DISPERSION
NRE-24048	SINGLE WALLED CARBON NANOTUBES(WATER) DISPERSION
NRE-24049	SWCNT WATER DISPERSION (>95%, DIAMETER: 1-2 nm, LENGTH: 5-15µm)
NRE-24050	SWCNTs DOPED WITH Ag NANOPOWDER IN WATER DISPERSION, 1WT%
NRE-24051	SWCNTs ETHANOL DISPERSION
NRE-24052	SWCNTs XYLENE DISPERSION
NRE-24053	SWCNTs-COOH WATER DISPERSION
NRE-24054	SWCNTs-COOH XYLENE DISPERSION
NRE-24055	SWCNTs-OH ETHANOL DISPERSION
NRE-24056	SWCNTs-OH WATER DISPERSION
NRE-24057	SWCNTs-OH XYLENE DISPERSION






CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-25001	CESIUM SELENIDE NANOPARTICLES	Cs ₂ Se, 99.9%, APS <100 nm
NRE-25002	COPPER INDIUM GALLIUM SELENIDE NANOPARTICLES	CuInGaSe ₂ , 99.9%, APS <100 nm
NRE-25003	GADOLINIUM METAL NANOPARTICLES	Gd, 99.9%, APS <100 nm
NRE-25004	GALLIUM ARSENIDE NANOPARTICLES	GaAs, 99.9%, APS <100 nm
NRE-25005	INDIUM ARSENIDE NANOPARTICLES	InAs, 99.9%, APS <100 nm
NRE-25006	INDIUM GALLIUM ARSENIDE NANOPARTICLES	InGaAs, 99.9%, APS <100 nm
NRE-25007	INDIUM SELENIDE NANOPARTICLES	In ₂ Se ₃ , 99.9%, APS <100 nm
NRE-25008	LANTHANUM III PHOSPHATE HYDRA NANOPARTICLES	LaPO ₄ • xH ₂ O, 99.9%, APS <100 nm
NRE-25009	POTASSIUM SELENIDE NANOPARTICLES	K ₂ Se, 99.9%, APS <100 nm
NRE-25010	RHODIUM METAL NANOPARTICLES	Rh, 99.9%, APS <100 nm
NRE-25011	RUBIDIUM SELENIDE NANOPARTICLES	Rb ₂ Se, 99.9%, APS <100 nm
NRE-25012	ZINC SELENIDE NANOPARTICLES	ZnSe, 99.9%, APS <100 nm
NRE-25013	NEODYMIUM FLUORIDE POWDER	NdF ₃ , 99.9%, APS <100 nm



Rare Earth Materials

Micro Powder


CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-26001	CESIUM SELENIDE POWDER	Cs ₂ Se, 99.9%, APS <40μm
NRE-26002	COPPER INDIUM GALLIUM SELENIDE POWDER	CuInGaSe ₂ , 99.9%, APS <40μm
NRE-26003	GADOLINIUM METAL POWDER	Gd, 99.9%, APS <40μm
NRE-26004	GALLIUM ARSENIDE POWDER	GaAs, 99.9%, APS <40μm
NRE-26005	INDIUM ARSENIDE POWDER	InAs, 99.9%, APS <40μm
NRE-26006	INDIUM GALLIUM ARSENIDE POWDER	InGaAs, 99.9%, APS <40μm
NRE-26007	INDIUM SELENIDE POWDER	In ₂ Se ₃ , 99.9%, APS <40μm
NRE-26008	LANTHANUM III PHOSPHATE HYDRA POWDER	LaPO ₄ • xH ₂ O, 99.9%, APS <40μm
NRE-26009	POTASSIUM SELENIDE POWDER	K ₂ Se, 99.9%, APS <40μm
NRE-26010	RHODIUM METAL POWDER	Rh, 99.9%, APS <40μm
NRE-26011	RUBIDIUM SELENIDE POWDER	Rb ₂ Se, 99.9%, APS <40μm
NRE-26012	ZINC SELENIDE POWDER	ZnSe, 99.9%, APS <40μm





CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-27001	ALUMINUM FOIL FOR LITHIUM ION BATTERY	Al, 99.9%
NRE-27002	ALUMINUM FOIL ROLL POLYMER BATTERY	Al, 99.9%
NRE-27003	ALUMINUM FOILS FOR ALUMINIUM STRIP	Al, 99.9%
NRE-27004	BATTERY PE SEPARATOR FOR LITHIUM ION BATTERY	99.90%
NRE-27005	CARBON CONDUCTIVE ADHESIVE TAPES	99.9%
NRE-27006	LITHIUM COBALT OXIDE NANOPARTICLES	LiCoO ₂ , 99.9%, <100 nm
NRE-27007	LITHIUM IRON PHOSPHATE	LiFePO ₄ , 99.9%
NRE-27008	LITHIUM MANGANESE OXIDE BATTERY MATERIAL	LiMn ₂ O ₄ , 99.9%
NRE-27009	LITHIUM NICKEL MANGANESE COBALT OXIDE	LiNiMnCoO ₂ , 99.9%
NRE-27010	LITHIUM TITANATE OXIDE POWDER	Li ₄ Ti ₅ O ₁₂ , >99%
NRE-27011	MESO CARBON MICROBEADS	99.9%
NRE-27012	NMP (N-METHYLPYRROLIDONE)	99.90%





CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-28001	ALUMINUM FOIL FOR LITHIUM ION BATTERY	Al, PURITY 99.9%
NRE-28002	ALUMINUM FOIL ROLL POLYMER BATTERY	Al, PURITY 99.9%
NRE-28003	ALUMINUM FOILS FOR ALUMINIUM STRIP	Al, PURITY 99.9%
NRE-28004	CARBON GRAPHITE PLATES	C, PURITY 99.9%
NRE-28005	COPPER FOIL	Cu, PURITY 99.9%
NRE-28006	COPPER FOIL ROLL	Cu, PURITY 99.9%
NRE-28007	COPPER FOIL SHEET ROLL	Cu, PURITY 99.9%
NRE-28008	COPPER FOIL SHEETS	Cu, PURITY 99.9%
NRE-28009	COPPER FOILS	Cu, PURITY 99.9%
NRE-28010	COPPER NICKEL FOIL	CuNi, PURITY 99.9%
NRE-28011	Cu FOIL	Cu, PURITY 99.9%
NRE-28012	Cu FOIL GRAPHENE	Cu, PURITY 99.9%
NRE-28013	GOLD FOIL	Au, PURITY 99.9%
NRE-28014	NICKEL FOIL	Ni, PURITY 99.9%
NRE-28015	PALLADIUM METAL FOIL	Pd, PURITY 99.9%
NRE-28016	TITANIUM METAL FOIL	Ti, PURITY 99.9%
NRE-28017	TUNGSTEN METAL FOIL	W, PURITY 99.9%



CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-58001	CHROMIUM ALUMINUM CARBIDE MAX PHASE POWDER	Cr ₂ AlC, PURITY 99.9%, APS<40μm
NRE-58002	MANGANESE ALUMINUM CARBIDE MAX PHASE POWDER	Mn ₂ AlC, PURITY 99.9%, APS<40μm
NRE-58003	MAX PHASE NIOBIUM ALUMINUM CARBIDE POWDER	Nb ₄ AlC ₃ , PURITY 99.9%, APS<40μm
NRE-58004	MAX PHASE TANTALUM ALUMINUM CARBIDE POWDER	Ta ₄ AlC ₃ , PURITY 99.9%, APS<40μm
NRE-58005	MAX PHASE TITANIUM ALUMINUM CARBIDE POWDER	Ti ₂ AlC, PURITY 99.9%, APS<40μm
NRE-58006	MAX PHASE TITANIUM SILICON CARBIDE POWDER	Ti ₃ SiC ₂ , PURITY 99.9%, APS<40μm
NRE-58007	MAX PHASE VANADIUM ALUMINUM CARBIDE POWDER	V ₄ AlC ₃ , PURITY 99.9%, APS<40μm
NRE-58008	MOLYBDENUM ALUMINUM CARBIDE MAX PHASE POWDER	Mo ₃ AlC ₂ , PURITY 99.9%, APS<40μm
NRE-58009	MOLYBDENUM GALLIUM CARBIDE MAX PHASE POWDER	Mo ₂ Ga ₂ C, PURITY 99.9%, APS<40μm
NRE-58010	MOLYBDENUM TITANIUM ALUMINUM CARBIDE MAX PHASE POWDER	Mo ₂ TiAlC ₂ , PURITY 99.9%, APS<40μm
NRE-58011	NIOBIUM ALUMINUM CARBIDE MAX PHASE POWDER	Nb ₂ AlC, PURITY 99.9%, APS<40μm
NRE-58012	TANTALUM ALUMINUM CARBIDE MAX PHASE POWDER	Ta ₂ AlC, PURITY 99.9%, APS<40μm
NRE-58013	TITANIUM ALUMINUM CARBIDE MAX PHASE POWDER	Ti ₃ AlC ₂ , PURITY 99.9%, APS<40μm
NRE-58014	TITANIUM ALUMINUM NITRIDE MAX PHASE POWDER	Ti ₂ AlN, PURITY 99.9%, APS<40μm
NRE-58015	TITANIUM TIN CARBIDE MAX PHASE POWDER	Ti ₂ SnC, PURITY 99.9%, APS<40μm
NRE-58016	VANADIUM ALUMINUM CARBIDE MAX PHASE POWDER	V ₂ AlC, PURITY 99.9%, APS<40μm



CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-59001	MXENE TITANIUM CARBIDE FILM	Ti2C, PURITY 99.9%, APS<40µm
NRE-59002	NIOBIUM CARBIDE MXENE POWDER	Nb2C, PURITY 99.9%, APS<40µm
NRE-59003	TITANIUM CARBIDE MXENE FILM (DELAMINATED)	Ti3C2, PURITY 99.9%, APS<40µm
NRE-59004	TITANIUM CARBIDE MXENE FILM (UNDELAMINATED)	Ti3C2, PURITY 99.9%, APS<40µm
NRE-59005	VANADIUM CARBIDE MXENE POWDER	V2C, PURITY 99.9%, APS<40µm



Targets & Wafers


Sputtering Target

CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-43001	ALUMINA SPUTTERING TARGET	Al ₂ O ₃ , PURITY: 99.99%
NRE-43002	ALUMINIUM CHROMIUM ALLOY SPUTTERING TARGET	AlCr, PURITY: 99.99%
NRE-43003	ALUMINIUM OXIDE SPUTTERING TARGET	Al ₂ O ₃ , PURITY: 99.99%
NRE-43004	ALUMINIUM OXINATE SPUTTERING TARGET	C ₂₇ H ₁₈ AlN ₃ O ₃ , PURITY: 99.99%
NRE-43005	ALUMINIUM SILICON SPUTTERING TARGET	AlSi, PURITY: 99.99%
NRE-43006	ALUMINIUM SPUTTERING TARGET	Al, PURITY: 99.99%
NRE-43007	ALUMINUM NITRIDE SPUTTERING TARGET	AlN, PURITY: 99.99%
NRE-43008	ANTIMONY SPUTTERING TARGET	Sb, PURITY: 99.99%
NRE-43009	BERYLLIUM SPUTTERING TARGET	Be, PURITY: 99.99%
NRE-43010	BISMUTH SPUTTERING TARGET	Bi, PURITY: 99.99%
NRE-43011	BORON CARBIDE SPUTTERING TARGET	B ₄ C, PURITY: 99.99%
NRE-43012	BORON NITRIDE SPUTTERING TARGET	BN, PURITY: 99.99%
NRE-43013	BORON SPUTTERING TARGET	B, PURITY: 99.99%
NRE-43014	CADMIUM SPUTTERING TARGETS	Cd, PURITY: 99.99%
NRE-43015	CALCIUM FLUORIDE SPUTTERING TARGET	CaF ₂ , PURITY: 99.99%
NRE-43016	CARBON SPUTTERING TARGET	C, PURITY: 99.99%
NRE-43017	CERIUM SPUTTERING TARGET	Ce, PURITY: 99.99%
NRE-43018	CHROMIUM SPUTTERING TARGET	Cr, PURITY: 99.99%
NRE-43019	COBALT CHROMIUM ALLOY SPUTTERING TARGET	Co ₃ Cr, PURITY: 99.99%
NRE-43020	COBALT IRON BORON SPUTTERING TARGET	CoFeB, PURITY: 99.99%
NRE-43021	COBALT IRON SPUTTERING TARGET	CoFe, PURITY: 99.99%
NRE-43022	COBALT OXIDE SPUTTERING TARGET	CoO, PURITY: 99.99%
NRE-43023	COBALT SPUTTERING TARGET	PURITY: 99.99%
NRE-43024	CONSTANTAN RESISTANCE ALLOY SPUTTERING TARGET	PURITY: 99.99%
NRE-43026	COPPER SPUTTERING TARGET	Cu, PURITY: 99.99%
NRE-43027	COPPER SULFIDE SPUTTERING TARGET	CuS, PURITY: 99.99%
NRE-43028	DYSPROSIUM SPUTTERING TARGET	Dy, PURITY: 99.99%
NRE-43029	ERBIUM SPUTTERING TARGET	Er, PURITY: 99.99%
NRE-43030	EUROPIUM SPUTTERING TARGET	Eu, PURITY: 99.99%
NRE-43031	GADOLINIUM SPUTTERING TARGET	Gd, PURITY: 99.99%
NRE-43032	GALLIUM NITRIDE SPUTTERING TARGET	GaN, PURITY: 99.99%
NRE-43033	GERMANIUM SPUTTERING TARGET	Ge, PURITY: 99.99%
NRE-43034	GOLD PALLADIUM SPUTTERING TARGET	AuPd, PURITY: 99.99%
NRE-43035	GOLD SPUTTERING TARGET	Au, PURITY: 99.99%
NRE-43036	GOLD TIN SPUTTERING TARGET	Au:Sn, PURITY: 99.99%
NRE-43037	HAFNIUM ALUMINIUM OXIDE SPUTTERING TARGET	HfAlO ₂ , PURITY: 99.99%
NRE-43038	HAFNIUM OXIDE SPUTTERING TARGET	HfO ₂ , PURITY: 99.99%
NRE-43039	INDIUM TIN OXIDE SPUTTERING TARGET	ITO, PURITY: 99.99%
NRE-43040	IRON SPUTTERING TARGET	Fe, PURITY: 99.99%
NRE-43041	LEAD SPUTTERING TARGET	Pb, PURITY: 99.99%
NRE-43042	LEAD ZIRCONATE TITANATE SPUTTERING TARGET	PZT, PURITY: 99.99%
NRE-43043	MAGNESIUM FLUORIDE SPUTTERING TARGET	MgF ₂ , PURITY: 99.99%
NRE-43044	MAGNESIUM FLUORIDE SPUTTERING TARGETS	MgF ₂ , PURITY: 99.9%
NRE-43045	MAGNESIUM SPUTTERING TARGET	Mg, PURITY: 99.9%
NRE-43046	MANGANESE OXIDE SPUTTERING TARGET	MnO, PURITY: 99.99%

Targets & Wafers



Sputtering Target

CODE NO.	PRODUCT NAME	DESCRIPTION
NRE-43047	MANGANESE SPUTTERING TARGET	Mn, PURITY: 99.99%
NRE-43048	MOLYBDENUM SPUTTERING TARGET	Mo, PURITY: 99.99%
NRE-43049	NICKEL SPUTTERING TARGET	Ni, PURITY: 99.99%
NRE-43051	PALLADIUM SPUTTERING TARGET	Pd, PURITY: 99.99%
NRE-43052	PLATINUM SPUTTERING TARGET	Pt, PURITY: 99.99%
NRE-43053	RHODIUM SPUTTERING TARGET	Rh, PURITY: 99.99%
NRE-43054	RUTHENIUM SPUTTERING TARGET	Ru, PURITY: 99.99%
NRE-43055	SILICON NITRIDE SPUTTERING TARGET	Si ₃ N ₄ , PURITY: 99.99%
NRE-43056	SILICON OXIDE SPUTTERING TARGET	SiO ₂ , PURITY: 99.99%
NRE-43057	SILICON SPUTTERING TARGET	Si, PURITY: 99.99%
NRE-43058	SILVER SPUTTERING TARGET	Ag, PURITY: 99.99%
NRE-43060	STAINLESS STEEL SPUTTERING TARGET	PURITY: 99.99%
NRE-43061	TIN SPUTTERING TARGET	Sn, PURITY: 99.99%
NRE-43062	TIN SULFIDE SPUTTERING TARGET	SnS, PURITY: 99.99%
NRE-43063	TITANIUM ALUMINIUM VANADIUM SPUTTERING TARGET	TiAlV, PURITY: 99.99%
NRE-43064	TITANIUM COPPER PALLADIUM ZIRCONIUM ALLOY SPUTTERING TARGET	TiCuPdZr, PURITY: 99.99%
NRE-43065	TITANIUM DIBORIDE SPUTTERING TARGET	TiB ₂ , PURITY: 99.99%
NRE-43066	TITANIUM OXIDE SPUTTERING TARGET	TiO ₂ , PURITY: 99.99%
NRE-43067	TITANIUM SPUTTERING TARGET	Ti, PURITY: 99.99%
NRE-43068	TRIPHENYLAMINE SPUTTERING TARGET	C ₅₇ H ₄₈ N ₄ , PURITY: 99.99%
NRE-43069	TUNGSTEN SPUTTERING TARGET	W, PURITY: 99.99%
NRE-43070	YTTERBIUM FLUORIDE SPUTTERING TARGET	YbF ₃ , PURITY: 99.99%
NRE-43071	YTTRIUM OXIDE SPUTTERING TARGET	Y ₂ O ₃ , PURITY: 99.99%
NRE-43072	YTTRIUM SPUTTERING TARGET	Y, PURITY: 99.99%
NRE-43073	ZINC OXIDE SPUTTERING TARGET	ZnO, PURITY: 99.5%
NRE-43074	ZINC SELENIDE SPUTTERING TARGET	ZnSe, PURITY: 99.99%
NRE-43075	ZINC SPUTTERING TARGET	Zn, PURITY: 99.99%
NRE-43076	ZINC SULFIDE SPUTTERING TARGET	ZnS, PURITY: 99.99%
NRE-43077	ZIRCONIUM COPPER ALUMINIUM SILVER ALLOY SPUTTERING TARGET	ZrCuAlAg, PURITY: 99.99%
NRE-43078	ZIRCONIUM SPUTTERING TARGET	Zr, PURITY: 99.99%







SILICON WAFER

Size (Inch)	2, 3, 4, 5, 6, 8
Diameter (mm)	50.8±0.3, 76.2±0.3, 100±0.5, 127±0.5, 154±0.5, 200±0.5
Growth method	CZ / FZ
Grade	PRIME / TEST / DUMMY GRADE
Thickness	180-1000 μm (+/-25) AS REQUIRED
Orientation	<100> / <110> / <111>
Type	N-TYPE/P-TYPE
Dopant	P-TYPE: BORON
Particle	N-TYPE: PHOSPHOROUS/ANTIMONY/ARSENIC
Oxygen	≤30 @0.3μm
Content	≤18 NEW PPMA
Carbon Content	≤1 NEW PPMA
Resistivity	1-10 ohm/cm
Others	TTV≤15μm, BOW≤65μm, WARP≤65μm
Surface	SINGLE SIDE POLISHED/DOUBLE SIDE POLISHED
Package	PACKED IN CASSETTE AND SEALED IN VACUUM BAG, 25 PCS / CASSETTE.



CODE NO.	PRODUCT NAME	DESCRIPTION
	SWCNT	
NRE-32005	HIGH PURITY SINGLE WALLED CARBON NANOTUBES	SWCNT, HIGH PURIFIED 98%
NRE-32009	SHORT HIGH PURITY SINGLE WALLED CARBON NANOTUBES	SWCNT, HIGH PURIFIED 98%
NRE-32013	INDUSTRIAL HIGH PURITY SINGLE WALLED CARBON NANOTUBES	SWCNT, PURITY >90%
	DWCNT	
NRE-33003	HIGH PURITY DOUBLE-WALLED CARBON NANOTUBES	DWCNTS, PURITY >98%
NRE-33004	DOUBLE-WALLED CARBON NANOTUBES	DWCNTS, PURITY >60%
NRE-33007	SHORT LENGTH DOUBLE-WALLED CARBON NANOTUBES	DWCNTS, PURITY >60%
	MWCNT	
NRE-34009	PURIFIED MWCNTs	OD:20-30 nm, LENGTH 10-20 µm >99WT%
NRE-34010	SHORT LENGTH MWCNTs	OD:20-30 nm, LENGTH 1-2 µm >99WT%
NRE-34014	INDUSTRIAL MWCNTs	OD:20-30 nm, LENGTH 10-20 µm >95WT%
	FUNCTIONALIZED	
NRE-34001	NH ₂ FUNCTIONALIZATION MWCNTs	OD:20-30 nm, LENGTH 10-20 µm >99WT%
NRE-34006	COOH FUNCTIONALIZED MWCNTs	OD:20-30 nm, LENGTH 10-20 µm >99WT%
NRE-34017	OH FUNCTIONALIZED MWCNTs	OD:20-30 nm, LENGTH 10-20 µm >99WT%
	GRAPHENE & GRAPHENE OXIDE	
NRE-39026	GRAPHENE POWDER	
NRE-39012	GRAPHENE OXIDE POWDER	
NRE-39009	GRAPHENE NANOPATELETS	
NRE-39021	REDUCED GRAPHENE OXIDE POWDER	
	FULLERENES	
NRE-41002	FULLERENE C ₆₀	PURITY >99%
NRE-41007	FULLERENE C ₇₀	PURITY >99%
NRE-41008	POLYHYDROXYLATED FULLERENE C ₆₀	PURITY >99%







Head Office
Plot No. 1010/121 , Sarita Vihar,
New Delhi-110076 ,India
Work Office
Village -Dhanora Jattan,
Kurukshetra Haryana-136132
Phone +91-786-402-0002
Phone +91-904-131-6296
E-mail - info@nanorh.com



Alberston , NY- 11507
USA
Phone +1-518-889-9730



Bridge St. , Glasgow
UK
Phone +44-141-846-0344

sales@nanorh.com

www.nanorh.com